INSTRUCTIONAL INFORMATION

INTENDED USE

Streck Cell Preservative is a solution that is formulated to preserve the white blood cells in peripheral blood samples without reducing the activity of antigenic sites. Samples treated with Streck Cell Preservative can be maintained for 7 days prior to processing and analysis by flow cytometric methods. Data has been developed supporting the use of Streck Cell Preservative for preserving the antigenic sites of lymphocytes in peripheral blood, lymph nodes, bone marrow and cord blood from normal donors and patient samples.

SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes in a blood sample. There are often situations where the sample cannot be analyzed immediately after being drawn, and must be preserved for future analysis. Streck Cell Preservative can preserve white blood cell integrity as well as the antigenic sites until processing and analysis can be performed. Blood samples can be transported to a distant site for analysis or stored until a larger batch of samples is accumulated for analysis.

Subsets of leukocytes can be distinguished on the basis of cell surface antigens using fluorescent antibodies and flow cytometry. Qualitative and quantitative changes in leukocyte subsets are used to identify and monitor immunodeficiency and hematologic diseases. Streck Cell Preservative is designed to preserve peripheral blood samples so that the qualitative and quantitative leukocyte subset characteristics do not differ from those of the original sample.

REAGENTS

Streck Cell Preservative is a cell preservative in a liquid medium.

PRECAUTIONS

For In Vitro Diagnostic Use. Studies have been conducted indicating that under specific conditions Cyto-Chex reagent will inactivate HIV. However, all human material must be handled with Universal Precaution protocols.

MSDS can be obtained on the homepage of www.streck.com, or by calling the following phone number, 1-800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

When stored at 18° to 30°C, Streck Cell Preservative is stable until the expiration date listed on the vial. It should not be frozen.

INDICATIONS OF PRODUCT DETERIORATION

1. Cloudiness or precipitate visible.
2. Patient samples do not show light scatter stability.
3. Patient samples do not react with fluorescent tagged antibodies.

INSTRUCTIONS FOR USE

The purpose of this diluent is to stabilize human blood cells up to seven days for subsequent processing and analysis by flow cytometry.

1. Add appropriate amount of patient sample to Streck Cell Preservative, equal parts reagent to sample, a 1:1 ratio (example: 1 ml patient sample to vial containing 1 ml Streck Cell Preservative).
2. Mix sample by inverting vial 3 times.
3.
 - a. For Becton-Dickinson Systems: The sample vial should be maintained at 2° to 10°C until use. This applies to storage and transportation.
 - b. For Beckman Coulter® Systems: The sample vial should be maintained at room temperature 18° to 30°C until use. This applies to storage and transport.
4. Mix sample vial thoroughly by hand inversion, at least 25 times.
5. Aliquot appropriate volume of sample into tubes for analysis.
6. Incubate with monoclonal antibody according to manufacturer's directions.
7. Process sample with red blood cell lyse/fix reagents according to manufacturer's directions.

PROCEDURES

1. **Instrument Procedure.** Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. **Fluorescent Antibody Procedure.** Use fluorescent antibodies according to manufacturer's instructions for patient and control samples.
3. **RBC Lysing Procedure.** Follow manufacturer's lysing reagent instructions.

LIMITATIONS

1. Do not dilute or add other components to Streck Cell Preservative.
2. Absolute counts must be adjusted for dilution factor.

EXPECTED RESULTS

Immunophenotyping performed on peripheral blood samples stored in Streck Cell Preservative should provide the same results as when performed on fresh specimens. NOTE: Absolute counts must be adjusted for the dilution factor.


REFERENCES

1. Turpen PF, Collins M. A reagent for stabilizing blood samples. *Am Clin Lab.* 1996 Sep; 15(8):30-1.
2. Schumacher MJ, Burkhead T. Stability of fresh and preserved fetal and adult lymphocyte cell surface markers. *J Clin Lab Anal.* 2000;14(6):320-6.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

Streck Cell Preservative was formerly labeled as Cyto-Chex®.

The brand and product names of the instruments are trademarks of their respective holders.

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