## INSTRUCTIONS FOR USE

### INTENDED USE

Cyto-Chex BCT is intended for collection and storage of blood specimens for immunophenotyping of white blood cells by flow cytometry. Recovery of lymphocyte subset cell markers of the HIV panel can be accomplished over a 14-day period following collection.

### SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes in a blood sample. Logistics, such as transport of samples from collection site to analysis location, create situations in which sample analysis cannot be performed immediately after collection. Sample preservation for future analysis provides a solution which addresses such situations.

Cyto-Chex BCT is a direct-draw blood collection tube containing an anticoagulant and preservative, which preserves the cell surface antigens of white blood cells (leukocytes) until processing and analysis can be performed.

Subsets of leukocytes can be distinguished on the basis of cell surface antigens using fluorescent antibodies and flow cytometry. Qualitative and quantitative changes in leukocyte subsets are used to identify and monitor immunodeficiency and hematologic diseases.

Cyto-Chex BCT is designed to preserve peripheral blood samples' qualitative and quantitative leukocyte subset characteristics.

### REAGENTS

Cyto-Chex BCT contains the anticoagulant, K<sub>3</sub>EDTA, and a cell preservative in a liquid medium.

### PRECAUTIONS

- For In Vitro Diagnostic Use.
- Do not freeze specimens collected in Cyto-Chex BCT as breakage could result.
- Do not use tubes after expiration date.
- Do not use tubes for collection of materials to be injected into patients.
- Product is intended for use as supplied. Do not dilute or add other components to Cyto-Chex BCT.
- Overfilling or under-filling of tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.
- CAUTION**
  - Glass has the potential for breakage; precautionary measures should be taken during handling.
  - All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection. Dispose of in accordance with federal, state and local regulations. Avoid contact with skin and mucous membranes.
  - Product should be disposed with infectious medical waste.
  - Remove and reinsert stopper by either gently rocking the stopper from side to side or by grasping with a simultaneous twisting and pulling action. A "thumb roll" procedure for stopper removal is not recommended, as tube breakage and injury may result.
- MSDS can be obtained on the homepage of [www.streck.com](http://www.streck.com), or by calling the following phone number, 1-800-843-0912, or by calling your local supplier.

### Prevention of Backflow

Since Cyto-Chex BCT contains chemical additives, it is important to avoid possible backflow from the tube. To guard against backflow, observe the following precautions:

- Keep patient's arm in the downward position during the collection procedure.
- Hold the tube with the stopper uppermost.
- Release tourniquet once blood starts to flow in the tube, or within 2 minutes of application.
- Tube contents should not touch stopper or the end of the needle during the collection procedure.

### STORAGE AND STABILITY

- When stored at 18 to 30°C, unused Cyto-Chex BCT is stable through expiration date.
- Do not freeze unfilled Cyto-Chex BCT. Proper insulation may be required for shipment during extreme temperature conditions.

### HIV Panel of Markers

- The HIV panel of markers including CD3, 4, 8, 16/56, 19 and 45 are stable in blood samples stored in Cyto-Chex BCT for up to 14 days at 18-22°C and up to 7 days at 30°C.
- CD3, 4, and 8 are stable for up to 48 hours at 37°C.
- CD3 and CD4 are stable for up to 72 hours at 37°C.

### Other CD Markers

Test results indicate CD markers 2, 5, 7, 10, 13, 20, 22, 23, 33, 34, 38, and HLA-DR are stable in blood samples stored in Cyto-Chex BCT for up to 7 days at 18-22°C.


### INDICATIONS OF PRODUCT DETERIORATION

- Cloudiness or precipitate visible.
- Patient samples do not exhibit typical light scatter patterns.
- Patient samples do not react with fluorescently-labeled antibodies.

## INSTRUCTIONS FOR USE

### Direct Draw Method

- Collect specimen by venipuncture according to CLSI H3-A6.
- Follow recommendations for order of draw outlined in CLSI H3-A6.
- Fill tube completely.
- Remove tube from adapter and immediately mix by gentle inversion 8 to 10 times. Inadequate or delayed mixing may result in inaccurate test results.

## PROCEDURES

- Cyto-Chex BCT should be at room temperature (18-25°C) and properly labeled for patient identification.
- Collect blood into the tube according to CLSI document H3-A6.
- After collection, transport and store tubes within the recommended temperature range to the laboratory for analysis.
- Perform immunophenotyping by flow cytometry in accordance with manufacturer's instructions.
 

Note:

  - Light scatter positions of cells stabilized with Cyto-Chex BCT may differ slightly from those of untreated cells.
  - Cyto-Chex BCT does not dilute blood samples; therefore, no dilution factor correction is necessary to obtain absolute count values.
  - As in the case with most clinical laboratory specimens, hemolysis, icteris and lipemia may affect the results obtained on blood samples preserved with Cyto-Chex BCT. In accordance with CAP Checklist for Flow Cytometry, the specimen condition should be noted on the laboratory report since the results may be suspect. Grossly hemolyzed samples should be rejected.

## LIMITATIONS

- Do not expose to temperatures less than 0°C or greater than 50°C.
- Samples drawn in other anticoagulants or preservatives may cause coagulation in Cyto-Chex BCT.

## EXPECTED RESULTS

Immunophenotypic analysis performed on peripheral blood samples stored in Cyto-Chex BCT will provide the same results as when performed on fresh specimens.

### CD4 + Cell Count Recovery

Healthy Donors	EDTA tube	Cyto-Chex BCT	
	Initial Cells/μl	Day 7 Cells/μl	Day 14 Cells/μl
1	597	684	641
2	1275	1421	1267
3	559	567	563
4	646	591	572
5	843	844	915
HIV+ Patients	Initial Cells/μl	Day 7 Cells/μl	Day 14 Cells/μl
1	709	760	756
2	511	561	565
3	285	299	300
4	327	340	373
5	83	73	89

## REFERENCES

- CLSI document H3-A6, Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.
- CDC Guidelines for Performing Single-Platform Absolute CD4+ T-Cell Determinations with CD45 Gating for Persons Infected with Human Immunodeficiency Virus, MMWR 52(RR-2), 2003.
- CLSI document H42-A2, Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline.
- CLSI document H43-A2, Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline.
- College of American Pathologists, Commission on Laboratory Accreditation, Flow Cytometry Checklist, 2003.

## ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

## GLOSSARY OF HARMONIZED SYMBOLS

EC REP	Authorized Representative in the European Community	REF	Catalog Number		Use By
LOT	Batch Code		Manufacturer		Temperature Limitation
IVD	In Vitro Diagnostic Medical Device		Consult Instructions For Use		Do not re-use
	Biological Risk	Glossary of symbols may contain symbols not used in the labeling of this product.			

U.S. Patents 5,196,182; 5,260,048; 5,460,797; 5,459,073; 5,849,517; 5,811,099; others pending

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## LITTERATUR

1. CLSI document H3-A6, Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.
2. CDC Guidelines for Performing Single-Platform Absolute CD4+ T-Cell Determinations with CD45 Gating for Persons Infected with Human Immunodeficiency Virus, MMWR 52(RR-2), 2003.
3. CLSI document H42-A2, Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline.
4. CLSI document H43-A2, Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline.
5. College of American Pathologists, Commission on laboratory Accreditation, Flow Cytometry Checklist, 2003.

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### RÉFÉRENCES

1. CLSI document H3-A6, Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.
2. CDC Guidelines for Performing Single-Platform Absolute CD4+ T-Cell Determinations with CD45 Gating for Persons Infected with Human Immunodeficiency Virus, MMWR 52(RR-2), 2003.
3. CLSI document H42-A2, Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline.
4. CLSI document H43-A2, Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline.
5. College of American Pathologists, Commission on Laboratory Accreditation, Flow Cytometry Checklist, 2003.

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1. CLSI document H3-A6, Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.
2. CDC Guidelines for Performing Single-Platform Absolute CD4+ T-Cell Determinations with CD45 Gating for Persons Infected with Human Immunodeficiency Virus, MMWR 52(RR-2), 2003.
3. CLSI document H42-A2, Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; pautas aprobadas.
4. CLSI document H43-A2, Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; pautas aprobadas.
5. College of American Pathologists, Commission on Laboratory Accreditation, Flow Cytometry Checklist, 2003.

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