## INSTRUCTIONS FOR USE

### INTENDED USE

A1c-Cellular is a bi-level, whole blood based, assayed control material for evaluating the accuracy and precision of HbA1c procedures.

### SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated and semi-automated procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of the instrument. A1c-Cellular's intact RBC formulation allows users to verify all steps of analysis procedures, including lysing of the RBC.

### REAGENTS

A1c-Cellular contains stabilized human red blood cells in a preservative medium.

### PRECAUTIONS

- For In Vitro Diagnostic Use.
- All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
- Control products are not to be used as calibrators.

### STORAGE AND STABILITY

A1c-Cellular should be stored at refrigerated temperatures (2° to 10° C) before and after opening. Under these conditions, unopened vials of A1c-Cellular are stable through the expiration date. After opening, A1c-Cellular is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2° to 10° C.

### INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration; however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

- Review control product package insert and the operating procedure of the instruments.
- Check expiration date of the product on the vial. Discard outdated products.
- Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Technical Services at 800-843-0912 or online at [www.streck.com](http://www.streck.com).

### INSTRUCTIONS FOR USE

- Remove vials of A1c-Cellular from the refrigerator and warm to room temperature (18° C to 30° C) for 15 minutes before use.
- To mix: **(Do not mix mechanically)**
  - Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
  - Mix by rapid inversion to ensure the cells are suspended.
  - Vials stored for an extended period of time may require extra mixing.
  - Gently invert the vials 8 to 10 times immediately before sampling.
- Refer to the instrument manual for the system in use for analyzing control materials.
- After sampling, return to refrigeration for maximum open-vial stability.

Barcodes on the cap-pierceable vials are intended for use with Bio-Rad® Variant II™ and Tosoh instruments. To use the barcode, align the barcode on the vial(s) with the barcode reader on the instrument. The barcode will label the control with "LC" for Level 1 and "HC" for Level 2. The Bio-Rad Variant instrument can use these barcodes to transfer the data into a QC log file. In order to not use the barcode function, face the barcode away from the instrument's barcode reader.

### LIMITATIONS

Do not use the product past the expiration date or if there is evidence of gross hemolysis or microbial contamination.

### EXPECTED RANGES

The values listed in this insert are derived from replicate analysis. The results listed on this insert are lot specific and were performed by the control manufacturer and/or independent laboratories using the same lot. The mean values recovered on your system should be within the stated expected range. Good laboratory practice includes establishment of a laboratory mean and standard deviation by replicate analyses. Assay ranges are provided only as a guide.

Values and ranges for instruments not listed on the Assay sheet must be established by the user. It is recommended that at least 10 consecutive analyses be performed on a properly calibrated instrument for each level to establish the "assay" mean.

### QUALITY CONTROL

Streck offers STATs®, an interlaboratory quality control program, and STATs-Link®, which provides internet access to STATs reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the STATs Department at 800-898-9563, or by fax at 402-333-7874.

### ORDERING INFORMATION

Please call our Customer Service Department toll free 800-843-0912 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

## GLOSSARY OF HARMONIZED SYMBOLS

ECIREP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	

Glossary of symbols may contain symbols not used in the labeling of this product.

The brand and product names of the instruments are trademarks of their respective holders.

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