INSTRUCTIONS FOR USE

INTENDED USE

ESR-Chex is an assayed bi-level control for evaluating the accuracy and precision of automated and manual Erythrocyte Sedimentation Rate (ESR) methods.

SUMMARY

Clinical laboratories are required to establish quality control programs for automated and semi-automated procedures. Daily use of ESR-Chex provides assurance that the erythrocyte sedimentation rate procedure is being performed properly.

PRINCIPLES

ESR-Chex is designed to be handled in the same manner as a patient specimen. There are several factors concerning the environment and techniques that may affect the outcome of erythrocyte sedimentation rate methods. When these are properly controlled, the values recovered with ESR-Chex will be within the Expected Ranges printed on the assay. Accurate analysis of patient specimens is assured.

REAGENTS

ESR-Chex contains stabilized human red blood cells in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.

STORAGE AND STABILITY

ESR-Chex is stable through the expiration date when stored at 2-10° C. After opening, ESR-Chex is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2-30° C.

INDICATIONS OF PRODUCT DETERIORATION

Discoloration of the product may be caused by overheating or freezing during shipping or storage. Gross hemolysis (darkly colored supernatant) may be indicative of product deterioration. However, moderately colored supernatant is normal and should not be confused with deterioration of the product.

INSTRUCTIONS FOR USE

1. INSTRUMENT PROCEDURE. Follow the instrument manufacturer's instructions for ESR testing for automated methods. Use the "Westergren 1 hr" setting for the Ves-Matic™ and Mini-Ves™ instruments.
2. MANUAL PROCEDURE. For manual testing, use ESR-Chex as you would a patient sample.
3. CAUTION. DO NOT remove sodium citrate or sodium chloride from tubes before using this control.

HANDLING INSTRUCTIONS

1. Remove vials from refrigerator and allow them to equilibrate to room temperature (20-30 minutes).
2. Vortex unopened vials for up to 60 seconds prior to first time use. Unopened vials stored at 2-10° C may require extra mixing. Mix previously opened vials through inversion and by vigorously rolling upright between palms until red cells are completely suspended. Continue to mix for 90 seconds.
3. Allow bubbles to disperse and re-mix by inversion prior to sampling.
4. Draw the sample immediately after thorough mixing is completed.
If mixed vials sit for more than 1 minute before drawing the sample, the vial must be remixed by repeating step 2. Incomplete mixing can invalidate both the sample drawn and the remaining product in the vial.
5. Follow the manufacturer's directions for filling the sedimentation rate tube for both automated and manual systems.
6. Wipe threads of vial and cap with clean tissue before closing. Recap the vial tightly.
7. We recommend storing open vials at room temperature (18-30° C). It will be easier to resuspend cells in vials stored at room temperature. Opened ESR-Chex vials are stable throughout the open-vial dating, as indicated on the assay sheet.
8. Previously opened vials should not require additional vortexing. Re-mix by inversion prior to sampling.
9. It is good laboratory practice to prepare a fresh sample of ESR-Chex daily to control the ESR procedure.

LIMITATIONS

ESR-Chex is to be used for Erythrocyte Sedimentation Rate testing only. It is not intended for controlling other hematology procedures.

EXPECTED RESULTS

The assay values are derived from replicate analysis on both automated and manual methods. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits. However, the control means established by the laboratory should fall within the Expected Range specified for the control.

If established means do not fall within the Expected Range, please evaluate all mechanical and physical factors that could affect the outcome, such as temperature, vibration, tube position and product expiration. Inadequate mixing of the product can also affect the outcome of the test, see Handling Instructions. Sample must be taken within 60 seconds of thoroughly mixing the vial. Assay an unopened vial of ESR-Chex. If the values are still outside the Expected Range, contact Streck Technical Services, 800-843-0912, or consult Technical Services at www.streck.com.


QUALITY CONTROL PROGRAM

Streck offers STATs®, an interlaboratory quality control program, and STATs-Link®, which provides internet access to STATs reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the STATs Department at 800-898-9563, or by fax at 402-333-7874. Additional information can be found online at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

The brand and product names of the instruments are trademarks of their respective holders.

U.S. Patents 5,895,760; 5,863,799; 5,888,822, 6,017,764, 6,051,433; 6,124,089; 6,159,682; 6,265,148; 6,342,391; 6,331,435; 6,531,321

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