## INSTRUCTIONS FOR USE

### INTENDED USE

ESR-Vacuum Tubes are evacuated tubes used for the collection of venous blood. ESR-Vacuum Tubes are used to transport and process blood for testing the Erythrocyte Sedimentation Rate (ESR) of whole blood in the clinical laboratory.

ESR analysis can be performed over a 72-hour period following collection and storage at 2-10° C in an ESR-Vacuum Tube.

### SUMMARY AND PRINCIPLES

The ESR test is a common laboratory test used to indicate patient conditions such as acute tissue damage, chronic infection and chronic inflammation. The ESR test measures the rate that erythrocytes settle in human plasma over a specified time period. ESR results are reported in millimeters that the plasma-erythrocyte interface falls in a designated perpendicular tube of anticoagulated whole blood per hour. ESR results become elevated due to the formation of erythrocyte rouleaux, which is caused mainly by alterations in plasma and erythrocyte factors.

There may be situations where the sample cannot be analyzed immediately after it is drawn and must be transported or preserved for future analysis. ESR-Vacuum Tubes preserve the integrity of the patient sample for sed rate analysis until processing and testing can be performed. Samples can be maintained from the time of blood collection for up to 72 hours prior to analysis when transported and/or stored at 2-10° C, or up to 4 hours at 18-30° C.

Safety Coated ESR-Vacuum Tubes\* are designed with an outer Mylar® coating that contains glass and blood specimens in the event of breakage, reducing the risk of injury and potential exposure to blood borne pathogens. Reports from an independent packaging consultant demonstrate that Safety Coated ESR-Vacuum Tubes are 10 times more impact resistant than glass collection tubes, significantly reducing accidental breakage.

ESR-Vacuum Tubes are available in 1.2ml and 2.0ml vacuum draw configurations, and are designed specifically for use with the ESR-100, ESR-8 and ESR-Auto Plus® automated sedimentation rate analyzers. Blood samples are collected in ESR-Vacuum Tubes and mixed with tri-sodium citrate to avoid coagulation. After the tube is mixed, it is placed in the ESR instrument. Each instrument automatically prints or downloads the ESR result at 30 minutes or 1 hour, depending on the program selected. 1.2ml ESR-Vacuum Tubes are also compatible with the ESR-10 Manual Rack, which produces a manual ESR result in 30 minutes.

### REAGENTS

ESR-Vacuum Tubes contain 3.2% buffered tri-sodium citrate (0.109M) additive.

### PRECAUTIONS

- For In Vitro Diagnostic Use.
- All patient blood samples should be treated as if they are capable of transmitting infections and should be handled with appropriate precautions. Avoid contact with skin and mucous membranes.
- Used tubes containing blood should not be disposed of in general waste, but should be disposed with infectious medical waste.
- The product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied, prior to direct draw or transfer procedures, invalidates the intended use of the product.
- Storage of tubes at or below 0° C may result in breakage.

### STORAGE AND STABILITY

ESR-Vacuum Tubes are stable through the expiration date when stored at 2° to 30° C. Samples can be maintained from the time of blood collection for up to 72 hours prior to analysis when transported or stored at 2-10° C, or up to 4 hours at 18-30° C.


### INDICATION OF PRODUCT DETERIORATION

- Inability to obtain expected ESR results may indicate product deterioration.
- Discoloration of the product may be caused by overheating or freezing during shipping or storage.
- Repeat the ESR run to confirm erroneous results; if the same result is obtained, contact Streck Technical Service at 800-843-0912 or online at [www.streck.com](http://www.streck.com).

### INSTRUCTIONS FOR USE

- Direct Draw Method
  - Follow the CLSI Approved Standard, "Procedure for the Collection of Diagnostic Blood Specimens by Venipuncture".
  - Insert the ESR-Vacuum Tube into the plastic holder and hold the tube so the cap is flush against the needle holder as the sample is being collected.
  - Angle the tube so the blood stream hits the glass wall before mixing with the citrate solution to minimize the formation of blood foam.
  - The vacuum will automatically draw the appropriate amount of blood into the tube. Watch for an air bubble to rise in the sample as the tube fills to indicate the draw is complete and remove the tube immediately.
  - The ideal fill level and acceptable fill range are indicated on the tube.
- Transfer Method
  - Blood samples drawn into an evacuated EDTA tube may be transferred into the ESR-Vacuum Tube. Do not remove anticoagulant from tube before transfer. Refer to CLSI guidelines for EDTA tube stability.
  - Mix the EDTA tube thoroughly by inverting the tube 6 to 8 times.
  - Remove the lids of the EDTA and ESR-Vacuum tubes. Transfer the sample from the EDTA tube into the ESR-Vacuum Tube up to the fill line.
  - The ideal fill level and acceptable fill range are indicated on the tube.
- Mix the citrate solution with the blood immediately after drawing or transferring by inverting the tube 6 to 8 times. The air bubble in the tube must reach the opposite end between every inversion. Hold the tube at a 35° angle to speed mixing.
- Additional Identification Labels
  - Additional identification labels placed on the tube must be placed outside the scanning area of the tube. Place any additional labels as close to the cap as possible. Leave a space of at least 3mm between the blood level and any additional labels.
  - ESR Tube Labeling Sleeves (Streck part number 240363) are available to provide an easy method of identifying samples without obstructing the blood column during analysis. Refer to ESR Tube Labeling Sleeve IFU for complete instructions.
    - Remove all 5 labeling sleeves from the sheet by folding along the perforated edges. Separate the

- individual labeling sleeves by folding along the perforated lines between each sleeve.
- Collect or transfer a valid blood specimen into an ESR-Vacuum Tube.
- Bend both sides of the labeling sleeve by creasing the tabs along the scored lines as shown in figure 1.
- Insert the tube through the hole in the labeling sleeve with the tabs pointing up around the rubber stopper as shown in figure 2.
- Apply the patient-specific identification label across both tabs of the labeling sleeve as shown in figure 3.

### LIMITATIONS

- For Single Use Only.
- The quantity of blood drawn may vary with altitude, temperature, barometric or venous pressure, tube age and technique. ESR-Vacuum Tubes are intended to be used with a standard vacutainer blood collection device. The use of a butterfly or other alternate blood collection device may affect the volume of blood drawn.
- Partial draw tubes (1.2ml) may fill more slowly than full draw tubes (2.0ml) of the same size. Under filling tubes will affect the final dilution of the sample, which may introduce variability into final results. Allow evacuated tubes to fill to completion.
- The level of vacuum within the 1.2ml High Altitude ESR-Vacuum Tube is slightly increased relative to the standard 1.2ml product. Use of this High Altitude tube at lower elevations (approximately 2500 ft. above sea level and below) may result in tubes filling well above the recommended product fill line.

ESR-Vacuum Tubes	1.2ml	2.0ml	1.2ml High Altitude	1.2ml Safety Coated
Nominal Capacity	1.3ml	2.2ml	1.6ml	1.3ml
Ideal Fill Level	60mm ± 5mm	100mm ± 8mm	60mm ± 5mm	60mm ± 5mm

### EXPECTED RESULTS

The ESR-Vacuum Tube should be filled with blood to the appropriate fill level indicated in the Limitations section above. When the ESR-Vacuum Tube is filled properly and tested in an ESR instrument according to the manufacturer's instructions, the results obtained on fresh blood samples will correlate with the modified Westergren method.


### REFERENCES

- Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Tubes and Additives for Venous Blood Specimen Collection; Approved Standard. CLSI Document H01.
- Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Reference and Selected Procedure for the Erythrocyte Sedimentation Rate (ESR) Test; Approved Standard. CLSI Document H02.
- Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard. CLSI Document H03.

### ORDERING INFORMATION

Please call our Customer Service Department toll free at 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

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### REFERENCES

1. Clinical and Laboratory Standards Institute (CLSI), ex-NCCLS : Tubes and Additives for Venous Blood Specimen Collection; Approved Standard. CLSI Document H01.
2. Clinical and Laboratory Standards Institute (CLSI), ex-NCCLS : Reference and Selected Procedure for the Erythrocyte Sedimentation Rate (ESR) Test; Approved Standard. CLSI Document H02.
3. Clinical and Laboratory Standards Institute (CLSI), ex-NCCLS : Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard. CLSI Document H03.

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### BIBLIOGRAFIA

1. Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Tubes and Additives for Venous Blood Specimen Collection; Approved Standard. CLSI Document H01.
2. Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Reference and Selected Procedure for the Erythrocyte Sedimentation Rate (ESR) Test; Approved Standard. CLSI Document H02.
3. Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard. CLSI Document H03.

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### NUORODOS

1. Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Tubes and Additives for Venous Blood Specimen Collection; Approved Standard. CLSI Document H01.
2. Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Reference and Selected Procedure for the Erythrocyte Sedimentation Rate (ESR) Test; Approved Standard. CLSI Document H01.
3. Clinical and Laboratory Standards Institute (CLSI), formerly NCCLS: Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard. CLSI Document H01.

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### REFERANSER

1. Clinical and Laboratory Standards Institute (CLSI), tidligere NCCLS: Tubes and Additives for Venous Blood Specimen Collection, Approved Standard. CLSI Document HO1.
2. Clinical and Laboratory Standards Institute (CLSI), tidligere NCCLS: Reference and Selected Procedure for the Erythrocyte Sedimentation Rate (ESR) Test; Approved Standard. CLSI Document HO2.
3. Clinical and Laboratory Standards Institute (CLSI), tidligere NCCLS: Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approval Standard. CLSI Document HO3.

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1. Clinical and Laboratory Standards Institute (CLSI) (antiguamente: NCCLS): Tubes and Additives for Venous Blood Specimen Collection; Approved Standard. Documento H01 de CLSI.
2. Clinical and Laboratory Standards Institute (CLSI) (antiguamente: NCCLS): Reference and Selected Procedure for the Erythrocyte Sedimentation Rate (ESR) Test; Approved Standard. Documento H02 de CLSI.
3. Clinical and Laboratory Standards Institute (CLSI) (antiguamente: NCCLS): Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard. Documento H03 de CLSI.

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