## INSTRUCTIONS FOR USE

### INTENDED USE

Para 12 Plus is an assayed hematology control for evaluating the accuracy and precision of hematology instruments that provide a white blood cell differential.

### SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation or manual protocols.

Use of stabilized cell preparations for controlling hematology instrumentation is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, Para 12 Plus will provide values within the expected range indicated on the assay sheet.

### REAGENTS

This whole blood reagent may contain any or all of the following: stabilized human or mammalian red blood cells, human, mammalian or simulated white blood cells and a platelet component in a preservative medium.

### PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. Control products are not to be used as calibrators.

### STORAGE AND STABILITY

Para 12 Plus is stable through the expiration date when stored at 2° to 10°C. After opening, Para 12 Plus is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2° to 10°C.

### INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of Para 12 Plus. Discard outdated products.
3. Assay an unopened vial of Para 12 Plus. If the values are still outside the Expected Range, contact Technical Services at 800-843-0912 or online at [www.streck.com](http://www.streck.com).

### INSTRUCTIONS FOR USE

#### Important Notes:

- ▶ The vial label barcode is only for use on the CELL-DYN Ruby.
  - ▶ CELL-DYN Enzymatic Cleaner Concentrate solution must be thoroughly removed from the analyzer Wash Block and/or Probe before running control products.
  - ▶ CELL-DYN 1700 instrument users: The CELL-DYN 1700 software will not recognize the leading zero digits in the lot number. You do not need to make any adjustments to the lot number when submitting your data to STATS.
1. Remove vials of control from the refrigerator and warm to room temperature (18° to 30°C) for 15 minutes before use.
  2. To mix: **(Do not mix mechanically)**
    - a. Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
    - b. Mix by rapid inversion to ensure the cells are suspended.
    - c. Vials stored for an extended period of time may require extra mixing.
    - d. Gently invert the vials 8 to 10 times immediately before sampling.
  3. Refer to the instrument manual for the system in use for analyzing control materials.
  4. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both vial and cap before replacing cap and returning to refrigeration.

### LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with Para 12 Plus. The white blood cell components simulate white blood cells in size, not morphology.

### EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. MCV and hematocrit values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Values and ranges for instruments not listed on the Assay sheet must be established by the user. It is recommended that at least 10 consecutive analyses be performed on a properly calibrated instrument for each level to establish the "assay" mean.

### QUALITY CONTROL

Streck offers *STATS*<sup>®</sup>, an interlaboratory quality control program, and *STATS-Link*<sup>™</sup>, which provides internet access to *STATS* reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the *STATS* Department at 800-898-9563, or by fax at 402-333-7874. Additional information can be found online at [www.streck.com](http://www.streck.com).

## ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

## GLOSSARY OF HARMONIZED SYMBOLS

ECREP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	

Glossary of symbols may contain symbols not used in the labeling of this product.

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