INSTRUCTIONS FOR USE

INTENDED USE

Retic-Chex II is an assayed control for evaluating the accuracy and precision of automated, semi-automated, and manual methods of reticulocyte counting. Retic-Chex Stain is used in conjunction with controls and patient samples for manual reticulocyte enumeration.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation or manual protocols.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument or manual methods, the whole blood control will provide values within the expected range indicated on the assay sheet.

REAGENTS

Retic-Chex II is composed of stabilized human red blood cells in a preservative medium. Retic-Chex Stain is filtered and buffered New Methylene Blue.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- Retic-Chex Stain - MSDS can be obtained on the homepage of www.streck.com, or by calling the following phone number, 1-800-843-0912, or by calling your local supplier.

= Irritant

www.streck.com

STORAGE AND STABILITY

Retic-Chex II is stable through the expiration date when stored at 2° to 10°C. After opening, Retic-Chex II is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2° to 10°C. DO NOT FREEZE.

Retic-Chex Stain can be shipped and stored at 2° to 30°C. After opening, Retic-Chex Stain is stable through the expiration date stated on the vial label.

INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the Expected Ranges:

- Review control product package insert and the operating procedure of the instrument.
- Check expiration date of product. Discard outdated products.
- Assay an unopened vial of product. If the values are still outside the Expected Range, call Technical Services at 800-843-0912 or online at www.streck.com.

INSTRUCTIONS FOR USE

Use Retic-Chex II Immediately After Removing From Refrigerator.

- Mix by gentle inversion between thumb and index finger until red blood cells are completely resuspended. Do not mix mechanically. Do not rub between palms of hands.
- Refer to appropriate procedure section below.
- Wipe threads of vial and cap with clean tissue before replacing cap. Recap vial.
- Return to refrigerator immediately.

MANUAL PROCEDURES

- Manual Method using a Miller Ocular:
 - Prepare dilution using equal drops of the control and Retic-Chex Stain.
 - Incubate at room temperature for a minimum of 15 minutes and no longer than 30 minutes.
 - Mix well. Prepare a smear and allow to dry.
 - Focus under 100x oil immersion optics on the feathered edge of the smear.
 - Count all reticulocytes contained in the large square area of the Miller Ocular. This includes reticulocytes in the smaller square.
 - Inspect and count several microscopic fields using the Miller Ocular. The total number of erythrocytes counted in the small squares should be at least 150 to assure reproducibility.
 - The percentage of reticulocytes is calculated according to the following formula:

$$\% \text{ reticulocytes} = \left(\frac{\text{Total reticulocytes in larger square}}{\text{Total RBCs in smaller square} \times 9} \right) \times 100$$
- Alternate Manual Method
 - Prepare the dilution and smears using procedure described in 1A, 1B and 1C.
 - Focusing under 100x oil immersion, count reticulocytes and mature red blood cells until 500 cells are accumulated.
 - Repeat procedure on another slide. The differences in reticulocyte counts in the two slide method must fall within the individual laboratory's established acceptable limits for duplicate analysis.
 - Divide the sum of reticulocytes counted on the two slides by ten to determine the reticulocyte percentage.

AUTOMATED PROCEDURES

- The user should follow the instrument manufacturer's instructions for performing automated reticulocyte counts.
- When using Retic-Chex II on the Abbott CELL-DYN® 3500/3700, it is recommended that analysis be performed 30 minutes after sample is added to the reagent.

EXPECTED RESULTS

When control values are within the Expected Range, patient values can be reported with assurance. The manual Retic-Chex II assay value is the average of replicate counts performed by using both dry methods. Automated and semi-automated assays are generated from properly calibrated instruments following instrument manufacturer's procedures. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the Expected Range specified for the control. Brilliant Cresyl Blue stain may be used with Retic-Chex II control, however, New Methylene Blue stain is the recommended stain according to CLSI Document H44-A2.

Values and expected ranges for instruments not listed on the Assay and Instructional Information sheet must

be established by the user. It is recommended that at least 10 consecutive analyses be performed on a properly calibrated instrument for each level to establish the "assay" mean.

QUALITY CONTROL PROGRAM

Streck offers STATS®, an interlaboratory quality control program, and STATS-Link™, which provides internet access to STATS reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the STATS Department at 800-898-9563, or by fax at 402-333-7874. Additional information can be found online at www.streck.com

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

ECREP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	

Glossary of symbols may contain symbols not used in the labeling of this product.

The brand and product names of the instruments are trademarks of their respective holders. U.S. Patent 5,432,089

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