INSTRUCTIONS FOR USE

INTENDED USE

Sickle-Chex is a positive and negative control for testing for the presence of sickling hemoglobin in solubility tests.

Sickle-Chex is compatible with the following sickle cell solubility kits: Streck *SICKLEDEX*[®], Dade[®] Behring Sickle-Sol[™] Solubility Test, Pacific Hemostasis SickleScreen[®] Sickling Hemoglobin Screening Kit, Chemo Diagnostic System Sickle-STAT and Columbia Calibre[®] Sickle Cell Reagent.

SUMMARY AND PRINCIPLES

Clinical laboratories are required to establish quality control programs for automated, semi-automated and manual procedures used for in vitro diagnostic testing of whole blood specimens.

Sickle-Chex is designed to be handled in the same manner as a patient specimen. There are several factors concerning the environment and technique that may affect the outcome of the test. When these are properly controlled, Sickle-Chex will produce the expected positive or negative result.

REAGENTS

Sickle-Chex is composed of stabilized human red blood cells in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specific by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. The product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. The product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied, invalidates any diagnostic use of the product.

STORAGE AND STABILITY

Open-vial stability is 100 days. After opening, Sickle-Chex is stable throughout the open-vial dating when stored at 2°C to 10°C, not to exceed the expiration date stated on the product vial.


INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. It is normal to observe a gradual darkening in the appearance of Sickle-Chex over the life of the product. This darkening is not an indication of product deterioration. If the recovered test results are not as expected:

1. Review control product package insert and the reagent test kit instructions.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Repeat the run to confirm erroneous results, if the same result is obtained, contact Technical Services at 800-843-0912 or online at www.streck.com.


INSTRUCTIONS FOR USE

1. Remove the vial of the control from the refrigerator and warm to room temperature (18° – 30°C) for 15 minutes before use.
2. To mix: (**Do not mix mechanically**)
 - a. Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
 - b. Mix by rapid inversion to ensure the cells are suspended.
 - c. Vials stored for an extended period of time may require extra mixing.
 - d. Gently invert the vials 8 to 10 times immediately before sampling.
3. Refer to the reagent test kit instructions for analyzing control and patient materials.
4. **Test kits that require a 20 µl sample should use 1 drop of Sickle-Chex control. Test kits that require a 50 µl sample should use 3 drops of Sickle-Chex control. For accurate delivery volume, the control vial **MUST** be inverted and held vertically directly over the test tube.**
5. After sampling, return to refrigeration for maximum open-vial stability. Wipe the threads of both vial and cap before replacing cap and returning to refrigeration.

EXPECTED RESULTS

The results of the Sickle-Chex control are read macroscopically. Incubate the control for the proper incubation time as indicated by the reagent test kit. The control results can be determined by holding the vial in front of the ruled lines in the example below or by using the viewing device included with the solubility test kit.

The positive test for sickling hemoglobin will result in a cloudy, turbid suspension and the ruled lines will not be visible through the vial. A negative control solution will appear virtually transparent and the ruled lines will be visible through the vial. Sickle-Chex controls are designed to verify the activity of the reagent. Negative patient results may not clear as quickly or as completely as the control.

REFERENCES

CONTROL - = negative control

CONTROL + = positive control

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

The brand and product names of the instruments are trademarks of their respective holders.

GLOSSARY OF HARMONIZED SYMBOLS

ECIREP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	

Glossary of symbols may contain symbols not used in the labeling of this product.

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