INSTRUCTIONAL INFORMATION

INTENDED USE

nRBC-Chex for ADVIA® is an assayed whole blood control designed to evaluate the accuracy and precision of the Siemens Healthcare Diagnostics Inc. hematology analyzers in the measurement of the nucleated red blood cell parameter. Refer to product assay sheet.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, the whole blood control will provide values within the expected range indicated on the assay sheet.

REAGENTS

The control (product) is composed of a mixture of stabilized human and animal blood cells. They are suspended in a solution containing a buffer and preservatives.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates the intended use of the product.
5. Control products are not to be used as calibrators.

INSTRUCTIONS FOR USE

1. Remove a vial of the control from the refrigerator and warm to room temperature (18° to 30°C) for 15 minutes before use.
2. To mix: (Do not mix mechanically.)
 - a. Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
 - b. Mix by rapid inversion to ensure the cells are suspended.
 - c. Vials stored for an extended period of time may require extra mixing.
 - d. Gently invert the vials 8 to 10 times immediately before sampling.
3. Refer to the instrument manual for the system in use for analyzing control materials. **Run nRBC-Chex for ADVIA in Patient Mode.**
4. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both vial and cap before replacing cap and returning to refrigeration.

STORAGE AND STABILITY

nRBC-Chex for ADVIA is stable through the expiration date when stored at 2° to 10°C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2° to 10°C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration; however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of nRBC-Chex for ADVIA on the vial. Discard outdated products.
3. Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Technical Service at 800-843-0912 or online at www.streck.com.

LIMITATIONS

A manual analysis of nucleated red blood cells cannot be accomplished with nRBC-Chex for ADVIA.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique. Code-outs may occur occasionally, however this does not indicate a failure of instrument or product unless it is observed excessively.

QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, and *STATS-Link*®, which provides Internet access to *STATS* reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the *STATS* Department at 800-898-9563, or by fax at 402-333-7874. Additional information can be found online at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

The brand and product names of the instruments are trademarks of their respective holders.

GLOSSARY OF HARMONIZED SYMBOLS

Glossary of symbols may contain symbols not used in the labeling of this product.

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