INSTRUCTIONS FOR USE**INTENDED USE**

CD-Chex CD34 is a stabilized preparation of human blood to be used as a complete process control when evaluating CD34 positive cells. It is intended to be used with BD Biosciences ProCOUNT™ Progenitor Cell Enumeration Kit, Beckman Coulter® Stem-Kit™, and with systems using the ISHAGE protocol for enumeration of CD34 cells.

SUMMARY AND PRINCIPLES

CD34 enumeration by flow cytometry provides a rapid and accurate assessment of the frequency of CD34 positive progenitor cells in samples from bone marrow, cord blood or peripheral blood from patients treated with hematopoietic growth factors. The ability to quantitate CD34 cells is useful in hematopoietic transplantation due to the fact that progenitor cells, required for successful engraftment, are found in the CD34 positive population.

CD34 positive progenitor cells can be distinguished on the basis of light scatter properties in conjunction with surface antigens. CD-Chex CD34 is designed to represent a blood sample containing CD34 positive cells having characteristics similar to progenitor cells: low side scatter properties, CD34 expression and low expression of CD45 (compared to lymphocytes). When stained with monoclonal antibodies for CD34 positive cell enumeration, CD-Chex CD34 control will provide reference values for CD34 positive cells within the ranges on the assay sheet.

REAGENTS

CD-Chex CD34 contains human leukocytes and erythrocytes in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.
5. CD-Chex CD34 products should not be used as a calibrator.

STORAGE AND STABILITY

CD-Chex CD34 is stable through the expiration date when stored at 2°C to 10°C. After opening, CD-Chex CD34 is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2°C to 10°C.

INDICATIONS OF PRODUCT DETERIORATION

If CD-Chex CD34 values are not within the expected range on the assay sheet or granulocytes show a loss of forward scatter (FSC):

1. Review the operating procedure of the instrument.
2. Assay an unopened vial of CD-Chex CD34. If values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or online at www.streck.com.
3. Clumping of the cell suspension indicates instability or deterioration, in which case the product should not be used.

INSTRUCTIONS FOR USE

CD-Chex CD34 is designed to be used with any standard CD34 enumeration protocol and has been evaluated with the ISHAGE, ProCOUNT and Stem-Kit CD34 enumeration protocols.

1. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18°C to 30°C) for 15 minutes before use.
3. To Mix:
 - a. Hold vial horizontally between the palms of the hands and roll the vial back and forth for 20 to 30 seconds.
 - b. Mix by gentle inversion at least 8-10 times until all cells are thoroughly suspended.
 - c. Gently invert the vial 5 times immediately before sampling.
4. Aliquot appropriate volume of control reagent into each tube and return to refrigeration to ensure maximum open vial stability.
5. Add recommended monoclonal antibodies according to manufacturer's instruction to each tube and mix gently.
6. Incubate according to antibody manufacturer's instructions.
7. Add recommended amount of RBC lysing agent according to manufacturer's instructions.
8. Analyze by flow cytometry in the same manner as patient samples.

RECOMMENDATION

CD-Chex CD34 has been evaluated with a variety of anti-CD34 mAb clones of various epitope class specificity. The following antibodies are recommended for use with CD-Chex CD34: HPCA2 (BD Biosciences), BIRMA-K3 (DakoCytomation) and 581 (Beckman Coulter Immunotech) which recognize class III epitopes of CD34.

LIMITATIONS

1. CD-Chex CD34 is not recommended for use with antibodies that target the class I and class II epitopes of CD34.
2. CD-Chex CD34 is designed to be used with a RBC lysing agent and may not provide results within the assay range if analyzed without RBC lysis.

EXPECTED RESULTS

Assay values provided for CD-Chex CD34 are derived from duplicate analyses on calibrated flow cytometers using the ISHAGE protocol and the ProCOUNT CD34 enumeration kit. If recovered values fall outside of this range, evaluate the instrument status, operator technique, and rule out product deterioration. Do not run patient specimens until the problem is resolved.


Since CD-Chex CD34 is manufactured from human blood, it contains white blood cells that bear cell surface antigens representative of normal leukocytes. For phenotyping white blood cells for antigens not listed on the assay sheet, it is recommended that at least five consecutive analyses be performed on a properly aligned and calibrated flow cytometer for each antibody in order to establish an "assay" mean.

QUALITY CONTROL

Streck offers STATs®, an interlaboratory quality control program, and STATs-Link®, which provides internet access to STATs reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the STATs Department at 800-898-9563, or by fax at 402-333-7874. Additional information can be found online at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

The brand and product names of the instruments are trademarks of their respective holders. U.S. Patents 5,196,182; 5,260,048; 5,459,073; 5,460,797; 5,811,099; 5,849,517

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