## INSTRUCTIONS FOR USE

### INTENDED USE

CD-Chex Plus is intended to be used as a quality control material for evaluating monoclonal antibody binding by flow cytometry. CD-Chex Plus control cells possess surface antigens detectable with monoclonal antibodies. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry they provide a reference for normal peripheral blood leukocytes. CD-Chex Plus is designed for use on the Becton Dickinson and Beckman Coulter® flow cytometry systems.

### SUMMARY/PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes present in a whole blood sample. Quantitation of peripheral blood leukocyte subsets has been shown to be useful in the evaluation of immunodeficiency, autoimmune, and oncologic disorders. CD-Chex Plus is a whole blood positive procedure control with assayed values used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance. CD-Chex Plus control cells mimic whole blood samples by possessing surface antigens detectable with fluorescent monoclonal antibodies by flow cytometry.

Subsets of leukocytes can be distinguished on the basis of light scatter properties as well as the presence of cell surface antigens. CD-Chex Plus is designed to represent normal peripheral blood leukocytes where lymphocytes, monocytes, and granulocytes are distinguishable on the basis of light scatter properties.

### REAGENTS

CD-Chex Plus contains stabilized human blood in a preservative medium.

### PRECAUTIONS

1. For In Vitro Diagnostic Use Only.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.
5. CD-Chex Plus should not be used as a calibrator.

### STORAGE AND STABILITY

CD-Chex Plus cells are stable through the expiration date when stored at 2-10° C. After opening, CD-Chex Plus is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2-10° C. DO NOT FREEZE.

### INDICATIONS OF PRODUCT DETERIORATION

If CD-Chex Plus values are not within the expected range on the assay sheet or granulocytes show a loss of forward scatter (FSC):

1. Review the operating procedure of the instrument.
2. Assay an unopened vial of CD-Chex Plus. If values are still outside the expected range, contact Technical Services, Streck, 800-843-0912 or www.streck.com.
3. Clumping of the cell suspension indicates instability or deterioration, in which case the reagent should not be used.

### INSTRUCTIONS FOR USE

1. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18-30° C) for 15 minutes before use.
3. To Mix:
  - a. Hold horizontally between palms of the hands and roll the vial back and forth for 20 to 30 seconds.
  - b. Mix by gentle inversion at least 8-10 times until all cells are thoroughly suspended.
  - c. Gently invert the vial 5 times immediately before sampling.

*Note: Vials stored for an extended period of time may require extra mixing.*
4. Aliquot appropriate volume of control reagent into each tube and return to refrigeration to ensure maximum open vial stability.
5. Add recommended monoclonal antibodies according to manufacturer's instruction to each tube and mix gently.
6. Incubate according to antibody manufacturer's instructions.
7. Add recommended amount of RBC lysing agent according to manufacturer's instructions.
8. Analyze by flow cytometry in the same manner as patient samples.

### LIMITATIONS

1. For optimum results, CD2, CD4, CD10 and CD19 should be evaluated using mAb conjugated to fluorochromes other than FITC. Antibody conjugates containing PE, RPE, RD1, ECD, and APC have been tested with CD-Chex Plus and are recommended for evaluation of CD2, CD4, CD10 and CD19.
2. CD2 recovery may be decreased when analyzed with Beckman Coulter® reagents and instrumentation.

### EXPECTED RESULTS

Mean and expected range values for each antigen are derived from replicate analyses on properly aligned and standardized flow cytometers. Results obtained using lysing systems, antibodies, and/or flow cytometers that differ from those used by Streck for value assignment may not fall within the expected range shown on the assay sheet. It is recommended that the user establish their own mean and acceptable range for each antigen analyzed by performing replicate analyses. Proper gating of all cell populations is necessary to obtain accurate results.

### QUALITY CONTROL

Streck offers STATS®, an interlaboratory quality control program, and STATS-Link®, which provides internet access to STATS reports, to all qualifying customers at no charge. If you are interested in more information or would like to participate, contact the STATS Department at 800-898-9563, or by fax at 402-333-7874. Additional information can be found online at www.streck.com.

## ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

## GLOSSARY OF HARMONIZED SYMBOLS

ECREP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	

Glossary of symbols may contain symbols not used in the labeling of this product.

The brand and product names of the instruments are trademarks of their respective holders. U.S. Patents 5,196,182; 5,260,048; 5,460,797; 5,459,073; 5,849,517; 5,811,099

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