

Technical Note

Detection of plasmid-mediated *ampC* β -lactamases using the Philisa® *ampC* ID Kit for identification of antibiotic resistance in gram-negative pathogens

Christopher M. Connelly, PhD¹, Stacey Morrow, MS, MT(ASCP)², Laura R. Porter¹, and Nancy D. Hanson, PhD²

Overview

AmpC β -lactamases are clinically important cephalosporinases that are resistant to most β -lactam antibiotics. AmpC enzymes are chromosomally encoded in many bacterial species and can be inducible and overexpressed as a consequence of mutation. Overexpression can lead to resistance to most β -lactam antibiotics. The occurrence of transmissible plasmids with acquired genes for AmpC β -lactamases often result in increased β -lactamase production, compared to chromosomally-expressed *ampC* genes. Additionally, plasmid-mediated AmpC β -lactamases can appear in organisms lacking or having low-level expression of a chromosomal *ampC* gene. Resistance due to plasmid-mediated AmpC enzymes can be broad in spectrum and often hard to detect. As such, it is useful to identify and discriminate between plasmid-mediated and chromosomally expressed AmpC β -lactamases. The Philisa® *ampC* ID Kit is a PCR-based molecular test that allows for multiplex identification of clinical isolates from six plasmid-mediated *ampC* gene families: MOX, DHA, ACC, EBC, FOX and CMY (see Table 1 for expanded list). An endogenous internal control is also included to reduce false negatives; it targets a conserved region common in gram-negative bacteria. Agarose gel detection is used to resolve PCR products and compare clinical samples against the external DNA controls. The Philisa *ampC* ID Kit can detect both plasmid-mediated and chromosomal *ampC* genes if the genes are not from the same chromosomal origin.

Table 1
Genes Identified with respective *ampC* ID primer sets

Primer Set	Gene(s)
MOX	MOX1-4, MOX8, CMY1, CMY8-11
ACC	ACC1, ACC2
FOX	FOX1-2, FOX4-9
DHA	DHA1, DHA2
CIT (CMY)*	CMY2, 4, 6, 7, 14-16, 18, 22, 25-44, 49, 53-56, 59, 60-62
EBC	ACT1-2, 8, 13 MIR1-3, 6-8
IC**	16S rRNA

* CIT refers to CMY-2-like genes that have their origin from *Citrobacter freundii*. Referenced as CMY in text.

** Internal control sequences are designed to detect *E. coli*, *Klebsiella* spp. and *Salmonella* spp.

The Philisa *ampC* ID Kit has been validated by extensive testing using previously characterized clinical isolates with the Philisa® Thermal Cycler and PhilisaFAST™ DNA Polymerase. Total PCR run time for this kit is 15 minutes, including hold times; the Philisa *ampC* ID Kit can rapidly screen test samples for the indicated gene families associated with antibiotic resistance. Optimization and validation studies demonstrated the Philisa *ampC* ID Kit's reliability when using PhilisaFAST and the Philisa Thermal Cycler. However, the Philisa *ampC* ID Kit is compatible with most DNA polymerases and thermal cycler platforms.

Materials and Methods

Testing of the Philisa *ampC* ID Kit was done at Creighton University School of Medicine, Department of Medical Microbiology and Immunology, Center for Research in Anti-Infectives and Biotechnology (C.R.A.B.). Data was provided with permission by Dr. Nancy Hanson.

PCR Amplification:

The Philisa *ampC* ID Kit (Catalog No.: 250026), Philisa Thermal Cycler (Catalog No.: 250000), PhilisaFAST DNA Polymerase (Catalog No.: 250024), and associated reagents were obtained from Streck, Inc. (Omaha, NE USA).

PCR was carried out as per the manufacturer's instructions for the Philisa *ampC* ID Kit: Hot start of 98 °C for 30 seconds, followed by 30 cycles [98 °C for 5 seconds, 58 °C for 10 seconds, and 72 °C for 7 seconds], and a final extension of 72 °C for 10 seconds. Rapid thermal cycling was carried out on the Philisa and the Bio-Rad C1000 Touch™. PCR amplicons were stained with ethidium bromide and resolved on a 2.5% agarose gel. A Molecular ImagerH Gel Doc XR+ System with Image Lab Software was used for PCR band detection and imaging.

