## INSTRUCTIONS FOR USE

### INTENDED USE

HQ-Chex is a whole blood glucose and hemoglobin control designed to evaluate the accuracy and precision of the HemoCue, Inc. glucose and hemoglobin analyzers.

### SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of instrument procedures that measure whole blood parameters. Daily use provides quality control data for confirming the precision and accuracy of instrument operation.

Use of stabilized cell preparations for controlling hematology and chemistry instrumentation is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, HQ-Chex will provide values within the expected range indicated on the assay sheet.

### REAGENTS

HQ-Chex contains stabilized human red blood cells and glucose in a preservative medium.

### PRECAUTIONS

- For In Vitro Diagnostic Use.
- CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at [www.streck.com](http://www.streck.com) for specific FDA required blood tests.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
- Assayed control products are not to be used as calibrators.
- SDS can be obtained at [www.streck.com](http://www.streck.com), by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

HQ-Chex is stable through the expiration date when stored at 2 °C to 10 °C. After opening, HQ-Chex is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. **DO NOT FREEZE.**

### INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

- Review control product package insert and the operating procedure of the instrument.
- Check expiration date of HQ-Chex. Discard outdated products.
- Assay an unopened vial of HQ-Chex. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

### PROCEDURES

**INSTRUMENT PROCEDURE.** Assay according to manufacturer's instructions for patient samples. Refer to assay values and ranges for the HemoCue analyzers.

### INSTRUCTIONS FOR USE

- Remove vials of control from the refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
- To mix (**do NOT mix mechanically or vortex**):

**For a video demonstration, visit [www.streck.com/mixing](http://www.streck.com/mixing).**

- Hold the vial vertically and roll each vial between the palms of the hands for 20-30 seconds.

- Continue to mix by holding the vial by the ends between the thumb and finger and rapidly inverting the vial 20 times end-over-end using a turning motion of the wrist.

- Visually inspect the bottom of the vial to ensure all cells are resuspended and sample immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5-10 times prior to sampling.
  - Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period.
- HemoCue analyzers - Do not fill the cuvette directly from the vial. Dispense a drop of the control onto a hydrophobic surface, for example a plastic film. Apply the cuvette to the drop.
  - After sampling, return to refrigeration for maximum open-vial stability. Wipe the threads of both the vial and cap before replacing cap and returning to refrigeration.

### EXPECTED RESULTS

The mean assay values provided are derived from replicate analyses on properly calibrated analyzers. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means

established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to instrument calibration, maintenance and operator technique.

### QUALITY CONTROL PROGRAM

Streck offers *STATS*, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or [statsdata@streck.com](mailto:statsdata@streck.com). Additional information can be found at [www.streck.com](http://www.streck.com).

### ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

### Rx Only

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [www.streck.com](http://www.streck.com).

The brand and product names of the instruments are trademarks of their respective holders.

See [www.streck.com/patents](http://www.streck.com/patents) for patents that may be applicable to this product.

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