STRECK 600

#P-117 APHL 2019

Detection of Mobilized Colistin Resistance (mcr) Genes by Multiplex **Real-Time PCR: Improving Surveillance** of an Emerging Threat

E. Alao, M.P. Torres, S. Cossette, C. Connelly Streck Research and Development, La Vista, NE

ABSTRACT

Recent increases in colistin-resistant infections led the CDC to launch an urgent public health response for mcr surveillance. Colistin is a last resort antibiotic that is more **Background** utilized due to increases in carbapenem-resistant infections. Since *mcr*-1 was first reported, more *mcr* gene variants have been identified, but few screening tools have been developed to rapidly detect mcr-positive samples. To improve surveillance for mcr genes, we describe a multiplex real-time PCR assay that detects mcr gene families 1 through 5 in less than 45 minutes.

This study utilizes sequence-specific primers and probes for the real-time PCR-based detection of mobilized colistin resistance mcr variants. An internal control (IC), Materials/ targeting a conserved region in Gram-negative bacteria, is also included in the multiplex mix to discriminate false negatives samples. Positive DNA controls are included with methods the multiplex assay. Data was generated using the Bio-Rad CFX96 Touch™ Real-time PCR Detection System.

Results

Conclusions

Streck ARM-D Kit, MCR provides a rapid amplification and detection strategy to monitor plasmid-mediated colistin resistance genes. The data demonstrate a sensitive and specific assay, with no observed cross-reactivity with previously characterized clinical isolates from Gram-negative organisms. The results demonstrate this assay can serve as a screening tool for surveillance of mcr-mediated colistin resistance, thereby improving antimicrobial stewardship practices to minimize mcr gene dissemination into the community

 Global increase in carbapenemase-producing Enterobacteriaceae has resulted in increased use of colistin.

 Colistin is considered a last resort antibiotic for multi-drug resistant Gram-negative bacteria.


were detected in 100% of samples and within 20 PCR cycles.

- Mobilized colistin resistance gene mcr-1 was first reported in 2015 in Escherichia coli strain isolated from pigs and humans in China.
- During the last 4 years mcr genes have been detected in over 30 countries and other types of bacteria such as Klebsiella pneumoniae and Salmonella enterica.
- According to the CDC, in the U.S. alone, over 50 human isolates have been reported across 19 states.

Figure 1: Tracking *mcr* in the U.S.

The following map shows where the *mcr* gene has been reported in U.S. human and food animal sources as of Nov. 2, 2018. Map source from the CDC website: https://www.cdc.gov/drugresistance/biggest-threats/tracking/mcr.html

BACKGROUND



The *mcr* real-time PCR multiplex assay is optimized to amplify *mcr* families 1 through 5. Amplification of serial dilutions of target controls generated PCR efficiencies ≥ 99% and correlation coefficients \geq 0.998. The sensitivity and specificity of the assay was evaluated using 90 clinical isolates and determined to be 100%. Internal control DNA

The results in **Figure 2** show the amplification of serial dilutions of ARM-D Kit, MCR Control Mix targets. Standard curves are shown for each target control. PCR efficiencies were over 99% for all the target controls and correlation coefficients were over 0.998. The amplification of the Multiplex Control Mix is shown in Figure 3 and respective C_a values for target controls are shown in **Table 1**.




Figure 2

Real-time PCR amplification of serially-diluted Multiplex Control Mix of Streck ARM-D Kit, MCR. Standard curves show corresponding efficiencies and correlation coefficients for each target control, respectively.







Figure 3

Real-time PCR amplification of Streck ARM-D Kit, MCR Multiplex Control Mix



Streck ARM¹ Control V Store at






<i>mcr</i> -like Gene Family	Average C _q Clinical Isolates	Positive Isolates (n)	Negative Isolates (n)	Sensitivity	Specificity
<i>mcr</i> -1	18.42 ± 0.34	8	82	100%	100%
mcr-2	19.36 ± 0.16	1	89	100%	100%
mcr-3	19.65 ± 0.4	2	88	100%	100%
mcr-4	17.5 ± 0.16	1	89	100%	100%
<i>mcr</i> -5*	19.27 ± 0.3	1	89	100%	100%
IC	14.56 ± 1.95	NA	NA	NA	NA
* All the DNA samples evaluated were extract	ad from clinical isolates excent for mcr-5				

* All the DNA samples evaluated were extracted from clinical isolates except for mcr-5 A contrived sample in gram negative bacteria matrix was used to evaluate sensitivity for this variant.

RESULTS

Out of the 90 isolates evaluated, 12 expressed *mcr* gene families. Representative amplification data of *mcr*positive isolates with respect to positive controls is shown in **Figure 4**. The *mcr* positive isolates amplified within 20 cycles of the PCR run (Table 2). The clinical isolates were correctly identified by the mcr assay and the sensitivity and specificity were 100%.



mcr-3-4-5: E=99.0%, R²=0.999





Table 1: C_a Values Streck ARM-D Kit, MCR Control Mix

<i>mcr</i> -like gene family	Fluorophore	Average C _q value of positive pontrol (n=2)
<i>mcr</i> -1	FAM	15.10 ± 0.13
mcr-2	HEX	16.47 ± 0.05
<i>mcr</i> -3-4-5	TEX615	16.18 ± 0
IC	CY5	14.64 ± 0.02






Figure 4

Real-time PCR amplification of *mcr* gene targets using Streck ARM-D Kits, MCR in clinical isolates (n= 90). Representative amplification plots of *mcr*-positive clinical isolates and positive control for each gene are shown. The internal control (IC) was amplified for all strains tested.

Table 2: Streck ARM-D Kt, MCR-Clinical Isolate Testing







4. Centers for Disease Control and Prevention Tracking the mcr gene-Updated on Nov. 2, 2018. https://www.cdc.gov/drugresistance/biggestthreats/tracking/mcr.html

MATERIALS & METHODS

Table 3: Gene targets covered by the Streck ARM-D Kit, MCR

ne Family	Gene Targets
<i>mcr</i> -1	mcr-1.1, 1.2, 1.3, 1.4, 1.5, 1.6, 1.7, 1.8, 1.9, 1.11, 1.12, 1.13, 1.14, 1.15
mcr-2	<i>mcr</i> -2.1
mcr-3	<i>mcr</i> -3.1, 3.2, 3.3, 3.4, 3.5, 3.6, 3.7, 3.8, 3.9, 3.10, 3.11, 3.12, 3.13, 3.14, 3.15, 3.16, 3.18, 3.19, 3.20, 3.21, 3.22, 3.23, 3.24, 3.25
mcr-4	mcr-4.1, 4.2, 4.3, 4.4, 4.5, 4.6
<i>mcr</i> -5*	<i>mcr</i> -5.1, 5.2, 5.3
IC	16S rRNA

The genes covered in Streck ARM-D Kit, MCR were differentiated using target-specific hydrolysis probes, chemically linked to different fluorescent dyes. The gene variants covered by the assay are shown in Table 3. PCR master mix preparation and PCR cycling protocol are shown in **Tables 4** and **5**, respectively. The efficiency and correlation coefficients were determined for each target by the amplification of serial dilutions of the Control Mix. The sensitivity and specificity of the assay was evaluated with DNA purified from *mcr*-negative and *mcr*-positive overnight bacterial cultures (n=90). Positive isolates for *mcr* gene targets were obtained from the CDC & FDA Antibiotic Resistance (AR) isolate bank panel with new or novel antibiotic resistance (AR Bank # 0346, 0349, 0350, 0493, 0494, 0495, 0496, 0497, 0538, 0539, 05400, 0635). DNA was isolated using the Qiagen® DNeasy® Blood and Tissue Kit or Exiprep Dx Bacteria Genomic DNA kit per the manufacturer's instructions. Duplicate reactions were run for all clinical isolates and controls.

Table 4: Master Mix Preparation			Table 5: PCR Cycling Conditions				
Source	Component	25 µL Reaction	Final Concentration	Step	Temp (°C)	Duration (sec)	Cycles
Lab supplied	Nuclease free water	9 μL	-	Enzyme Activation	98	30	Hold
Streck 2X Supermix	-	12.5 μL	1X	Denature	98	5	
reck ARM-D Kit	10X PCR Mix	2.5 μL	-	Anneal	60	10	30
Distribute Master Mix into PCR wells or tubes as appropriate before sample		Extension/Capture	72	10			
ab supplied or reck ARM-D Kit	Template Unknown or NTC -or- Template-Control	1μL	_				

CONCLUSIONS

The Streck ARM-D Kit, MCR:

- Identifies and differentiates mcr-1, mcr-2 and mcr-3-4-5 gene families with 100% sensitivity and 100% specificity.
- Can detect up to 48 out of 51 *mcr* allelic variants within the described gene families.
- Represents a surveillance tool to track the spread of mobilized colistin resistance.

REFERENCES & ACKNOWLEDGEMENTS

1. Sun J, Zhang H, Liu YH, Feng Y, 2018. Towards Understanding *mcr*-like Colistin Resistance. Trends Microbiol. 26(9): 794-808.

2. Caniaux I, van Belkum A, Zambardi G, Poirel L, Gros MF, 2017. mcr: Modern Colistin Resistance. Eur J Clin Microbiol Infect Dis. 36(3): 415-420.

3. Karaiskos I, Souli M, Galani I, Giamarellou H, 2017. Colistin: Still a lifesaver for the 21st century? Expert Opin Drug Metab Toxicol. 13(1): 59-71.

We would like to thank Streck for support of this work. We acknowledge and thank Jeremiah Athmer, Nancy Hanson (Creighton University) and Stacey Morrow (Creighton University) for their help with clinical isolates. We are also very grateful to Zach Friesen and Elizabeth Jorgenson for their assistance with the preparation of this poster.

