INSTRUCTIONS FOR USE

INTENDED USE

Streck Cell Preservative™ is a solution formulated to preserve the white blood cells in peripheral blood samples without reducing the activity of antigenic sites. Samples treated with Streck Cell Preservative can be maintained for up to 7 days prior to processing and analysis by flow cytometric methods.

SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes in a patient sample. There are often situations where the sample cannot be analyzed immediately after being obtained, and must be preserved for future analysis. Streck Cell Preservative can stabilize white blood cell integrity and the antigenic sites of leukocytes in peripheral blood, lymph nodes, bone marrow and cord blood until processing and analysis can be performed.^{1,2} Samples can be transported to a distant site or stored until a larger batch of samples is accumulated for analysis.

Subsets of leukocytes can be distinguished on the basis of cell surface antigens using fluorescent antibodies and flow cytometry. Qualitative and quantitative changes in leukocyte subsets are used to identify and monitor immunodeficiency and hematologic diseases. Streck Cell Preservative is designed to preserve the sample's qualitative and quantitative leukocyte subset characteristics for up to 7 days.

REAGENTS

Streck Cell Preservative is a general use reagent that contains the anticoagulant K₂EDTA and a cell preservative in a liquid medium.

PRECAUTIONS

1. To the best of our knowledge, unused product does not require any special disposal. However, each facility must determine proper disposal methods to comply with federal, state and local regulations.
2. Avoid contact with skin and mucous membranes.
3. Do not ingest.
4. Do not use reagent after expiration date.
5. Product is intended to be used as supplied. Do not dilute or add other components to Streck Cell Preservative.
6. Do not freeze specimens stored in glass vials, as breakage could result.
7. Upon exposure to biological specimens, this product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
8. SDS can be obtained at www.streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

1. Unopened, Streck Cell Preservative is stable until the date listed on the vial when stored at 2 °C to 30 °C. After initial vial opening, Streck Cell Preservative is stable for 30 days when stored at 2 °C to 30 °C.
2. Short-term storage at 2 °C to 40 °C is acceptable for unopened Streck Cell Preservative for up to 14 days.

INDICATIONS OF PRODUCT DETERIORATION

1. Cloudiness or precipitate visible.
2. If indications of product deterioration occur, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

The purpose of this reagent is to stabilize human blood cells up to 7 days for subsequent processing and analysis by flow cytometry.

Sample Dilution and Storage

1. Add appropriate amount of patient sample to Streck Cell Preservative, equal parts reagent to sample, a 1:1 ratio (example: 1ml patient sample to vial containing 1ml Streck Cell Preservative).
2. Mix sample by inverting vial 3 times.
3. a. For **BD Biosciences systems**: The sample vial should be maintained at 2 °C to 10 °C until use. This applies to storage and transport.
b. For **Beckman Coulter® systems**: The sample vial should be maintained at 18 °C to 30 °C until use. This applies to storage and transport.

Processing and Analysis

1. Mix sample vial thoroughly by hand inversion, at least 25 times.
2. Aliquot appropriate volume of sample into tubes for analysis.
3. Incubate with monoclonal antibodies according to manufacturer's directions.
4. Process sample with red blood cell lyse/fix reagents according to manufacturer's directions.

PROCEDURES

1. **Instrument Procedure.** Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. **Fluorescent Antibody Procedure.** Use fluorescent antibodies according to manufacturer's instructions.
3. **RBC Lysing Procedure.** Follow manufacturer's lysing reagent instructions.

LIMITATIONS

Absolute counts must be adjusted for dilution factor.

EXPECTED RESULTS

Immunophenotyping performed on peripheral blood samples stored in Streck Cell Preservative should provide equivalent results as when performed on fresh specimens.

REFERENCES

1. Turpen PB, Collins M. A reagent for stabilizing blood samples. *American Clinical Laboratory* 1996; 15(8):30-1.
2. Schumacher MJ, Burkhead T. Stability of fresh and preserved fetal and adult lymphocyte cell surface markers. *Journal of Clinical Laboratory Analysis* 2000; 14(6):320-6.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at www.streck.com.

Streck Cell Preservative was formerly labeled as Cyto-Chex®.

The brand and product names of the instruments are trademarks of their respective holders.

See www.streck.com/patents for patents that may be applicable to this product.

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