INSTRUCTIONS FOR USE

INTENDED USE

UA-Cellular[®] for IQ is an assayed cellular urine control for evaluating the accuracy and precision of automated procedures that measure urinary sediment parameters.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated procedures that measure urinary sediment parameters. Daily use of this urine cellular control provides quality control data for confirming the precision and accuracy of instrument operation.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, the urine cellular control will provide values within the expected range indicated on the assay sheet.

REAGENTS

UA-Cellular for IQ contains the following: stabilized human red blood cells and white blood cells and other inert particles in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. **CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at www.streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates the intended use of the product.
5. Control products are not to be used as calibrators.
6. SDS can be obtained at www.streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

UA-Cellular for IQ is stable through the expiration date when stored at 2 °C to 10 °C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the bottle. Discard outdated products.
3. Assay an unopened bottle of the control. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

1. Remove a bottle of the control from the refrigerator and warm at room temperature (18 °C to 30 °C) for 15 minutes before use.

Note: Remove foil liner from new bottles of control.

2. Refer to the instrument manual for the system specific instructions for analyzing patient samples. Do not utilize the control files of the instrument for UA-Cellular for IQ.

3. **For a video demonstration, visit www.streck.com/mixing.**

Mix by inverting end-over-end 30 times over approximately 30 seconds, using quick flicks of the wrist. Invert the container and verify that all material is off the bottom of the bottle. Do not mix mechanically. Due to the nature of the particles being analyzed, the control settles rapidly. Sample tubes should be prepared and analyzed immediately after mixing.

4. Place sample tubes for Level 1 in position 1 of the sample rack and Level 2 in position 2 (no more than 2 control tubes should be tested in any single rack due to settling of the particles during the instrument sampling process). Initiate analysis immediately. Repeat steps 2-4 for any additional testing of the control.
5. Return the control bottles to refrigeration for maximum open-vial stability.

LIMITATIONS

A manual sediment analysis cannot be accomplished with this product. The components simulate urinary sediment in size, not morphology.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

QUALITY CONTROL PROGRAM

Streck offers *STATS*[®], an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at www.streck.com.

The brand and product names of the instruments are trademarks of their respective holders.

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