## INSTRUCTIONS FOR USE

### INTENDED USE

A1c-Cellular<sup>®</sup> Linearity is a multi-level whole-blood based, assayed control material for evaluating the accuracy and precision of HbA1c procedures over the instrument reportable range.

### SUMMARY AND PRINCIPLES

Regulatory agencies mandate that laboratories substantiate their test methods throughout the reportable range for patient test results. Circumstances which may call for verification include installation, major preventive maintenance, unusual trends in control performance, a complete change of reagents or whenever recommended by the instrument manufacturer. A1c-Cellular Linearity contains HbA1c% concentrations which span typical patient reportable ranges, allowing the user to comply with these guidelines.

Assay values are derived from replicate testing on instruments operated and maintained according to manufacturer's instructions.

### REAGENTS

A1c-Cellular Linearity contains stabilized human red blood cells in a preservative medium.

### PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. **CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied unless otherwise directed. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product, unless pretreatment is recommended by individual instrument manufacturers.
5. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

A1c-Cellular Linearity should be stored at refrigerated temperatures (2 °C-10 °C) before and after opening. Under these conditions, unopened vials of A1c-Cellular Linearity are stable through the expiration date. After opening, A1c-Cellular Linearity is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C-10 °C.

### INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration; however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review product package insert and the operating procedure of the instruments.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Remix and rerun. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

## INSTRUCTIONS FOR USE

### Preparation


1. Perform routine maintenance, quality control procedures, calibration and precision checks according to your laboratory protocol.
2. Remove vials from refrigerator and warm to room temperature (18 °C-30 °C) for 15 minutes before use.

### Mixing and Handling

1. To mix: (**Do NOT mix mechanically or vortex.**)

**For a video demonstration, visit [streck.com/mixing](http://streck.com/mixing).**

- a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
- d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).

**Important Note for Roche cobas 6000, Roche cobas c 311, Siemens Dimension Series, VITROS 4600 / 5600 / XT 7600 and VITROS 5,1 FS Users:**

**REFER TO THE ASSAY SHEET FOR REQUIRED DILUTION PREPARATION.**

## Analysis Protocol for A1c-Cellular Linearity

1. Analyze all levels from lowest to highest vial number in Patient Mode.
2. Run each vial 3 times, mixing before each sampling event.
3. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both the vial and cap before replacing cap and return to refrigeration.
4. Instruments with assay specific values should verify that the mean values recovered are within the expected ranges. Key data into the product-specific *STATS*<sup>®</sup> Data Input Form located at streck.com and email the completed form to statsdata@streck.com, or login to your *STATS-Link*<sup>®</sup> account to submit your data.

### LIMITATIONS

Do not use the product past the expiration date or if there is evidence of gross hemolysis or microbial contamination.

### EXPECTED RESULTS

Values recovered should fall within the suggested tolerance limits on the accompanying assay. However, lab to lab variation in instrument performance, instrument calibration and operator technique may yield differences in results obtained. It is good laboratory practice that an individual lab establish its own means and limits for each parameter.

1. If your mean values are outside of your laboratory's established acceptance limits, review maintenance, QC procedures and calibration checks according to your lab protocol and rerun A1c-Cellular Linearity.
2. If the mean values remain outside your established limits, contact your instrument service representative. **DO NOT ADJUST YOUR INSTRUMENT.**
3. **Some levels may exceed the lower and/or upper limits of the reportable range provided by the manufacturer. Refer to the manufacturer operator's manual for specifications and performance claims.** It is important to note that terminology varies; CLIA<sup>1,2</sup> refers to the analytical range as the "Reportable Range" while CLSI<sup>3</sup> standards use the term "Analytical Measuring Interval" (AMI) and CAP<sup>4</sup> uses the term "Analytical Measurement Range" (AMR). AMR and AMI correspond to the CLIA "Reportable Range."

### REFERENCES

1. Centers for Medicare & Medicaid Services. Medicare, Medicaid, and CLIA Programs. Laboratory requirements relating to quality systems and certain personnel qualifications. Final Rule. Fed Regist 2003; 16: 3640-3714.
2. Centers for Medicare & Medicaid Services State Operations Manual, Appendix C - Survey procedures and interpretive guidelines for laboratories and laboratory services; Subpart K; Rev. 140, 2015.
3. Clinical and Laboratory Standards Institute, EP6-A, Evaluation of the linearity of quantitative measurement procedures: a statistical approach. Approved Guideline 2003.
4. College of American Pathologists, Commission on Laboratory Accreditation. Laboratory General Checklists. 2009.

### QUALITY CONTROL PROGRAM

Streck offers *STATS*, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

### ORDERING INFORMATION

Please call our Customer Service Department 800-843-0912 for assistance. Additional information can be found online at streck.com.

### Rx Only

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

Streck  
7002 S. 109 Street, La Vista, NE 68128 USA

MediMark<sup>®</sup> Europe  
11, rue Emile Zola, BP 2332  
38033 Grenoble Cedex 2, France

350596-12  
2020-03

**Barcodes on the cap-pierceable vials are intended for use with Bio-Rad D-10 and Tosoh instruments.**

To use the barcode, align the barcode on the vial(s) with the barcode reader on the instrument. The barcode will label the control with the appropriate level and lot number. In order to not use the barcode function, face the barcode away from the instrument's barcode reader.

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### QUELLENGABEN

- Centers for Medicare & Medicaid Services. Medicare, Medicaid, and CLIA Programs. Laboratory requirements relating to quality systems and certain personnel qualifications. Final Rule. Fed Regist 2003; 16: 3640-3714.
- Centers for Medicare & Medicaid Services State Operations Manual, Appendix C - Survey procedures and interpretive guidelines for laboratories and laboratory services; Subpart K; Rev. 140, 2015.
- Clinical and Laboratory Standards Institute, EP6-A, Evaluation of the linearity of quantitative measurement procedures: a statistical approach. Approved Guideline 2003.
- College of American Pathologists, Commission on Laboratory Accreditation. Laboratory General Checklists. 2009.

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