INSTRUCTIONS FOR USE

INTENDED USE

UA-Cellular[®] Complete is an assayed chemistry and cellular urine control for evaluating the accuracy and precision of the procedures that measure urinary sediment and chemistry parameters on the following instruments and testing methods: Sysmex UF-1000i Automated Urine Particle Analyzer; Siemens Clinitek Atlas Automated Urine Chemistry Analyzer utilizing the Clinitek Atlas 10 Reagent Pak; the Arkray AUTION HYBRID AU-4050 fully-automated integrated urine analyzer, utilizing the Uriflet S 9HA Test Strips; Siemens Clinitek Status line of automated chemistry strip reader analyzers; Manual Reading of Urine Refractive Index; Manual Reading of the Siemens Multistix 10SG Reagent Strips and Siemens Clinitest hCG Pregnancy Test.

The list of assayed parameters includes:

Sysmex UF-1000i: RBC (/μl), WBC (/μl), Epithelial (/μl), Cast, Bacteria (/μl), Crystals)

Arkray AU-4050: RBC (/μl), WBC (/μl), Epithelial (/μl), Cast, Bacteria (/μl), Crystals)

Arkray AU-4050 with Uriflet S, 9HA Test Strips: Glucose (mg/dL), Bilirubin (As Measured), Ketones (mg/dL), Specific Gravity (As Measured), Blood (As Measured), pH (As Measured), Protein (mg/dL), Urobilinogen (EU/dL), Nitrite (As Measured), Leukocytes (As measured), Color (As Measured), Turbidity (As Measured)

Siemens Clinitek Atlas with 10 Reagent Pak: Glucose (mg/dL), Bilirubin (As Measured), Ketones (mg/dL), Specific Gravity (As Measured), Blood (As Measured), pH (As Measured), Protein (mg/dL), Urobilinogen (EU/dL), Nitrite (As Measured), Leukocytes (As measured), Color (As Measured), Clarity (As Measured)

Siemens Clinitek Status System line with Multistix 10SG Reagent Strips: Glucose (mg/dL), Bilirubin (As Measured), Ketones (mg/dL), Specific Gravity (As Measured), Blood (As Measured), pH (As Measured), Protein (mg/dL), Urobilinogen (EU/dL), Nitrite (As Measured), Leukocytes (As measured)

Siemens Multistix 10SG Reagent Strips: Glucose (mg/dL), Bilirubin (As Measured), Ketones (mg/dL), Specific Gravity (As Measured), Blood (As Measured), pH (As Measured), Protein (mg/dL), Urobilinogen (EU/dL), Nitrite (As Measured), Leukocytes (As measured)

Siemens Clinitest hCG Pregnancy Test Cassette: Pregnancy hCG (As Measured)

Manual Confirmatory Testing: Refractive Index (nD)

SUMMARY

UA-Cellular Complete is designed specifically for use with: Sysmex UF-1000i Automated Urine Particle Analyzer; Siemens Clinitek Atlas Automated Urine Chemistry Analyzer utilizing the Clinitek Atlas 10 Reagent Pak; the Arkray AUTION HYBRID AU-4050 fully-automated integrated urine analyzer, utilizing the Uriflet S 9HA Test Strips; Siemens Clinitek Status line of automated chemistry strip reader analyzers; Manual Reading of Urine Refractive Index; Manual Reading of the Siemens Multistix 10SG Reagent Strips and Siemens Clinitest hCG Pregnancy Test.

The comprehensive tri-level UA-Cellular Complete control serves to align six stabilized urine particle types (i.e., red blood cells, white blood cells, casts, crystals, bacteria and epithelials) with the 10 urine chemistry parameters, plus hCG. The complexity of the UA-Cellular Complete is not duplicated by today's current commercial urine chemistry and urine cellular/particle type controls.

The overview below describes the tri-level configuration of UA-Cellular Complete:

Level PC Control – Serves as the positive chemistry control for the urine chemistry parameters that do not align with a corresponding particle type (i.e. glucose, bilirubin, ketones and urobilinogen).

Level PCM Control – Serves as the positive chemistry control for the urine chemistry parameters that align with a corresponding urine particle type as follows:

- Positive "Blood" is aligned with true human "RBC"
- Positive "Leukocyte Esterase" is aligned with true human "WBC"
- Positive "Nitrite" is aligned with true "Bacteria"
- Positive "Protein" is aligned with the WBC and RBC cellular components

Level NC Control – Serves as the normal control for the qualitative and semi-qualitative urine chemistry measurands. Protein is present in trace amounts. This level offers a lower range of RBCs and WBCs than Level PCM to comprehensively test the instruments' quantitative reporting range for these parameters. This level serves as the positive control for Siemens Clinitest hCG Pregnancy Test.

Laboratories require assayed material for quality control of automated and manual procedures that measure urinary sediment and chemistry parameters. Daily use of this urine control provides quality control data for confirming the precision and accuracy of instrument operation and lab technologist proficiency. UA-Cellular Complete is intended, at minimum, to be run by the clinical laboratory in conjunction with daily instrument and/or manual testing start-up to collectively verify the performance attributes (i.e., calibration, precision and accuracy) of the distinct instrument platforms and manual confirmatory tests, prior to running clinical patient specimens. Recovery of the assayed measurands within the stated expected ranges serve to confirm that the instrument platforms and manual tests are reporting correctly.

Use of stabilized cell preparations and chemistry analytes for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, UA-Cellular Complete will provide values within the expected range indicated on the assay sheet.

REAGENTS

UA-Cellular Complete contains the following: stabilized mammalian red blood cells and white blood cells, stabilized bacteria, and simulated urine sediments in a preservative medium. Analyte levels are adjusted with appropriate chemicals.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All human source material should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates the intended use of the product.
5. Control products are not to be used as calibrators.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

UA-Cellular Complete is stable through the expiration date when stored at 2 °C to 10 °C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. Do not Freeze this product. Store UA-Cellular Complete in its original, tightly closed bottle. Protect from light exposure. UA-Cellular Complete is light sensitive. Do not store in direct sunlight.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedures of the instrument.
2. Check expiration date of the product on the bottle. Discard outdated products.
3. Assay an unopened bottle of the control. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

1. This product should be treated the same as a patient specimen. Refer to the instrument manual for the system specific instructions or the product package insert for the manual confirmatory tests regarding analyzing patient samples.
2. Remove the control from the refrigerator and warm at room temperature (18 °C to 30 °C) for 15 minutes prior to use and protect from light exposure. **NOTE:** Remove foil liner from the new bottles of control.
3. **For a video demonstration, visit streck.com/mixing.** Mix by inverting end-over-end 30 times over approximately 30 seconds, using quick flicks of the wrist. Invert the container and verify that all material is off the bottom of the bottle. Do not mix mechanically. Due to the nature of the particles being analyzed, the control settles rapidly. Sample tubes should be prepared and analyzed immediately after mixing.

4. **Analyze the sample on the instrument within 15 seconds of re-suspending the bottle and transferring to a clean test tube. IMPORTANT:** Excessive foam or bubbles at the surface of the solution may cause sampling errors. If present, remove with a clean pipette.
5. Discard the used solution tubes after use.
6. Promptly return the bottle to refrigeration.

Sysmex UF-1000i additional instructions:

1. Perform an Auto-Rinse immediately prior to and after running the controls.
2. Refer to the instrument manual for the system specific instructions for analyzing patient samples. Do not utilize the control files of the instrument for UA-Cellular Complete.

LIMITATIONS

A manual sediment analysis cannot be accomplished with this product. Some of the components only simulate urinary sediment in size, not morphology.

UA-Cellular Complete must be well mixed prior to each use to ensure reproducible results. UA-Cellular Complete is specifically formulated for the instruments and manual tests listed in the INTENDED USE section.

UA-Cellular Complete is light sensitive and should be protected from light to achieve intended recovery.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Do not make any additions, chemical or otherwise, to UA-Cellular Complete. Altering the control could seriously affect results due to chemical and cellular compatibility and/or dilution.

QUALITY CONTROL PROGRAM

Streck offers **STATS[®]**, an interlaboratory quality control program, to all customers at no charge. For more information, contact the **STATS** Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

The brand and product names of the instruments are trademarks of their respective holders.

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