INSTRUCTIONS FOR USE

INTENDED USE

CD-Chex CD117[®] Plus is intended to be used as a quality control material for evaluating CD117, CD25 and CD71 monoclonal antibody binding by flow cytometry. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry, they provide a reference value for the abnormal cells found in certain types of hematopoietic neoplasms. CD-Chex CD117 Plus is designed for use on BD Biosciences and Beckman Coulter flow cytometry systems. **This product and the markers provided on the assay have not been cleared by the U.S. Food and Drug Administration for In Vitro Diagnostic use. This product and the values provided are For Research Use Only. Not for use in diagnostic procedures.**

SUMMARY AND PRINCIPLES

CD-Chex CD117 Plus is designed to represent abnormal peripheral blood leukocytes similar to a hematolymphoid neoplastic patient sample^{1,2,3}. CD-Chex CD117 Plus possesses surface CD117, CD25 and CD71 that are detectable with fluorescent monoclonal antibodies by flow cytometry. Abnormal leukocytes are distinguishable from normal leukocytes on the basis of light scatter properties and a low level of CD45 expression. CD-Chex CD117 Plus is a positive procedural assayed control used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance.

REAGENTS

CD-Chex CD117 Plus contains stabilized human blood and cells of human origin in a preservative medium.

PRECAUTIONS

1. CD-Chex CD117 Plus is For Research Use Only.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates the use of the product.
5. CD-Chex CD117 Plus should not be used as a calibrator.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

CD-Chex CD117 Plus is stable through the expiration date when stored at 2 °C to 10 °C. After opening, CD-Chex CD117 Plus is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. DO NOT FREEZE.

INDICATIONS OF PRODUCT DETERIORATION

If CD-Chex CD117 Plus values are not within the expected range on the assay sheet:

1. Review control product package insert, assay and operating procedures of the instrument.
2. Check expiration date of the product vial. Discard outdated products.
3. Clumping of the cell suspension indicates instability or deterioration, in which case the reagent should not be used.
4. Assay an unopened vial of CD-Chex CD117 Plus. If the values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

1. Follow instrument manufacturer's instructions for instrument compensation and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. Mixing Procedure (**mechanical mixing by vortex or rotator is not recommended**):
For a video demonstration, visit streck.com/mixing.
 - a. Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- b. Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.

- c. Aliquot immediately after mixing.
 - d. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.
- Note: Vials stored for an extended period of time may require extra mixing.
4. Return control reagent to refrigeration after sampling to ensure maximum open-vial stability.
 5. Add recommended monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.

Note: A negative staining control is recommended due to the heterogeneous expression of CD117, CD25 and CD71.

6. Incubate according to antibody manufacturer's instructions.
7. Add recommended amount of RBC lysing agent and follow manufacturer's instructions.
8. Analyze by flow cytometry by gating the abnormal cells and/or a total leukocyte "Live"[™] gate from a CD45/SSC or FSC/SSC plot. See Figure 1.

GATING

The most common gating strategy used in neoplastic cell assessment is to gate on the abnormal cells and then determine the CD marker percent positivity using a negative staining control¹⁻⁴. Abnormal cells are generally located with either a FSC/SSC plot or a CD45/SSC plot; although other gating strategies can be employed¹⁻⁴. A less frequently employed method is to gate all of the leukocytes in a "Live"[™] gate, then determine the CD marker percent positivity of the white blood cells.


FIGURE 1. Gating methods used to obtain CD117, CD25 and CD71 reference values.

In CD-Chex CD117 Plus, the abnormal cells (RED OVALS) are CD45^{dim}/FSC^{heterogeneous}/SSC^{intermediate}. These cells are positioned in the monocyte region of a FSC/SSC plot (A) or in the blast region of a CD45/SSC plot (B). In CD-Chex CD117 Plus, the "Live" gates (BLACK RECTANGLES) are the total leukocytes, which include all of the cells, normal and abnormal, excluding debris.

LIMITATIONS

1. For optimum results, CD117 should be evaluated using a mAb conjugated to fluorochromes other than PerCP-Cy5.5.
2. CD117 evaluated with mAb 95C3 when using lyse other than ammonium chloride will result in sub-optimal recovery.
3. CD25 recovery is reduced when using the Beckman Coulter mAb B1.49.9 conjugated to FITC.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on properly compensated flow cytometers⁴. The assay values are obtained using common flow cytometry reagents. See the assay for additional limitations or specific instructions for reagents.

Use of alternate gating strategies not specified in these instructions may result in values outside the published range. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own means and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different reagents, laboratory protocols, instrument calibration, maintenance, and operator technique. Data collected from interlaboratory quality control programs can be used as a cumulative approach when calculating ranges.

REFERENCES

1. Gorczyca W. Flow Cytometry in Neoplastic Hematology, Morphologic-Immunophenotypic Correlation-Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria S., Wood B. Flow Cytometry in Evaluation of Hematopoietic Neoplasms; A Case-Based Approach. CAP Press, Northfield, IL USA, 2012.
3. Nguyen D., Diamond L.W., Braylan R.C. Flow Cytometry in Hematopathology, A Visual Approach to Data Analysis and Interpretation-Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute H43-A2 Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline - Second Edition. Clinical and Laboratory Standards Institute, Wayne, PA USA, 2009.

ADDITIONAL RESOURCES

van Dongen J.J., Lhermitte L., Böttcher S., Almeida J., van der Velden V.H., Flores-Montero J., et al. EuroFlow Consortium. EuroFlow Antibody Panels for Standardized n-Dimensional Flow Cytometric Immunophenotyping of Normal, Reactive and Malignant Leukocytes. Leukemia. 2012 Sep; 26(9):1908-75.

Wood B.L., Arroz M., Barnett D., DiGiuseppe J., Greig B., Kussick S.J., et al. 2006 Bethesda International Consensus Recommendations on the Immunophenotypic Analysis of Hematolymphoid Neoplasia by Flow Cytometry: Optimal Reagents and Reporting for the Flow Cytometric Diagnosis of Hematopoietic Neoplasia. Cytometry B Clin Cytom. 2007; 72 Suppl 1:S14-22.

QUALITY CONTROL PROGRAM

Streck offers *STATS*[®], an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

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BIBLIOGRAFIA

1. Gorczyca W. Flow Cytometry in Neoplastic Hematology, Morphologic-Immunophenotypic Correlation-Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria S, Wood B. Flow Cytometry in Evaluation of Hematopoietic Neoplasms; A Case-Based Approach. CAP Press, Northfield, IL USA, 2012.
3. Nguyen D., Diamond L.W., Braylan R.C. Flow Cytometry in Hematopathology, A Visual Approach to Data Analysis and Interpretation-Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute H43-A2 Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline - Second Edition. Clinical and Laboratory Standards Institute, Wayne, PA USA, 2009.

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van Dongen J.J., Lhermitte L., Böttcher S., Almeida J., van der Velden V.H., Flores-Montero J., et al. EuroFlow Consortium. EuroFlow Antibody Panels for Standardized n-Dimensional Flow Cytometric Immunophenotyping of Normal, Reactive and Malignant Leukocytes. Leukemia. 2012 Sep; 26(9):1908-75.

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