## INSTRUCTIONS FOR USE

### INTENDED USE

Sperm-Chex® Post VC is a positive/negative control intended for monitoring post vasectomy sperm counts performed manually using commercially available counting chambers to validate the presence or absence of sperm in seminal fluids.

### SUMMARY AND PRINCIPLES

Microscopic evaluation of post vasectomy semen is a crucial part of the vasectomy process. The post vasectomy sperm count validates that the vasectomy procedure was successful. When handled like a patient sample, Sperm-Chex Post VC controls will provide a positive sperm count sample and a negative sperm count sample to validate the manual post vasectomy sperm count procedure.

### REAGENTS

Sperm-Chex Post VC control contains stabilized mammalian sperm in a preservative medium and a blank or negative control diluent.

### PRECAUTIONS

1. For in Vitro Diagnostic Use.
2. Product contains dilute formalin.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. SDS can be obtained at [streck.com](http://streck.com), by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

Sperm-Chex Post VC control is stable through the expiration date when stored at 2 °C to 10 °C. Cap the vial tightly and store upright. After opening, Sperm-Chex Post VC is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C.

### INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected results may indicate product deterioration.

To troubleshoot:

1. Review the control product package insert Instructions For Use and Limitations sections.
2. Check the expiration date of Sperm-Chex Post VC. Discard outdated product.
3. Assay an unopened vial of Sperm-Chex Post VC. If the expected results are not achieved, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

## INSTRUCTIONS FOR USE

For a video demonstration, visit [streck.com/mixing](http://streck.com/mixing).

1. Remove the controls from the refrigerator. It is not necessary to warm the controls to room temperature before using.
2. To mix:
  - a. Pipette any liquid from the cap and add to the contents of the vial before mixing. Do not invert vial during mixing.
  - b. **Vortex on medium speed for 3 to 5 pulses of 3 to 5 seconds each, verifying that a vortex is obtained. The pellet on the bottom of the vial should be dislodged and the clear liquid may turn turbid during the mixing process. Do not mix the product in any other manner. Other mixing methods will not properly re-suspend the sperm.**
3. Perform the post vasectomy count according to the hemocytometer package insert, Makler counting chamber or your established laboratory protocols.
4. After sampling, return the vial to refrigeration for maximum open-vial stability.

### LIMITATIONS

1. Analysis of Sperm-Chex Post VC cannot be performed on automated cell counting systems.
2. Sperm-Chex Post VC control is a positive or negative sperm count only and cannot be used for morphology, motility, or viability testing.

### EXPECTED RESULTS

The presence of sperm is identified as "present" for positive or "not present" for negative, does not have a recovery range and is determined visually using hemocytometers or Makler chambers on representative vials. Always perform quality control according to good laboratory practice, laboratory director's requirements, and any regulatory or accredited requirements.

### ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at [streck.com](http://streck.com).

### Rx Only

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [streck.com](http://streck.com).

See [streck.com/patents](http://streck.com/patents) for patents that may be applicable to this product.

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