Retic-Chex® II

Retic-Chex® Stain

CE

STRECK 

INSTRUCTIONS FOR USE

INTENDED USE

Retic-Chex® II is an assayed control for evaluating the accuracy and precision of automated, semi-automated, and manual methods of reticulocyte counting. Retic-Chex® Stain is used in conjunction with controls and patient samples for manual reticulocyte enumeration.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation or manual protocols.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument or manual methods, the whole blood control will provide values within the expected range indicated on the assay sheet.

REAGENTS

Retic-Chex II is composed of stabilized human red blood cells in a preservative medium. Retic-Chex Stain is filtered and buffered New Methylene Blue.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. **CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at www.streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. Retic-Chex Stain – SDS can be obtained at www.streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

Retic-Chex II is stable through the expiration date when stored at 2 °C to 10 °C. DO NOT FREEZE. After opening, Retic-Chex II is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C and Retic-Chex Stain is stable through the expiration date stated on the vial label. Retic-Chex Stain can be shipped and stored at 2 °C to 30 °C.

INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the Expected Range:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of product. Discard outdated products.
3. Assay an unopened vial of product. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

Use Retic-Chex II Immediately After Removing From Refrigerator.

1. To Mix: (Do NOT mix mechanically or vortex. Do not rub between palms of hands.)

For a video demonstration, visit www.streck.com/mixing.

- a. Mix by holding vial between thumb and index finger and gently inverting the vial 20 times end-over-end until the red blood cells are completely resuspended.

- b. Examine the bottom of the vial to ensure that the product is completely resuspended and repeat mixing as necessary.
- c. Mixing must be repeated upon removing the sample from the refrigerator for the entire open-vial time period.

2. Refer to appropriate procedure section below.

3. Wipe threads of vial and cap with clean tissue before replacing cap, and return to refrigeration immediately.

Note: Dropper tip may be carefully removed to pipette directly from vial. Carefully replace dropper tip prior to recapping vial.

MANUAL PROCEDURES

1. Manual Method using a Miller Ocular:

- a. Prepare dilution using equal drops of the control and stain.
- b. Incubate at room temperature for a minimum of 15 minutes and no longer than 30 minutes.
- c. Mix well. Prepare a smear and allow to dry.
- d. Focus under 100x oil immersion optics on the feathered edge of the smear.
 1. Count all reticulocytes contained in the large square area of the Miller Ocular. This includes reticulocytes in the smaller square.
 2. Inspect and count several microscopic fields using the Miller Ocular. The total number of erythrocytes counted in the small squares should be at least 150 to assure reproducibility.
 3. The percentage of reticulocytes is calculated according to the following formula:

$$\% \text{ reticulocytes} = \left(\frac{\text{Total reticulocytes in larger square}}{\text{Total RBCs in smaller square} \times 9} \right) \times 100$$

2. Alternate Manual Method

- a. Prepare the dilution and smears using procedure described in 1A, 1B and 1C.

b. Focus under 100x oil immersion optics where there is uniform distribution of the red cells.

1. Tally separately reticulocytes and mature red blood cells in successive fields until 500 total cells (number of reticulocytes and red cells) are accumulated.
2. Repeat the procedure on another slide. The differences in reticulocyte counts in the two-slide method must fall within the individual laboratory's established acceptable limits for duplicate analysis.
3. Divide the sum of reticulocytes counted on the two slides by 10 to determine the reticulocyte percentage.

AUTOMATED PROCEDURES

1. The user should follow the instrument manufacturer's instructions for performing automated reticulocyte counts.
2. When using Retic-Chex II on the Abbott CELL-DYN® 3500/3700, it is recommended that analysis be performed 30 minutes after sample is added to the reagent.

EXPECTED RESULTS

When control values are within the Expected Range, patient values can be reported with assurance. The manual Retic-Chex II assay value is the average of replicate counts performed by using both dry methods. Automated and semi-automated assays are generated from properly calibrated instruments following instrument manufacturer's procedures. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the Expected Range specified for the control.

Brilliant Cresyl Blue stain may be used with Retic-Chex II control, however, New Methylene Blue stain is the recommended stain according to CLSI Document H44-A2¹. Retic-Chex II is compatible with any brand of New Methylene Blue stain.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2², it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

REFERENCES

1. Clinical and Laboratory Standards Institute, H44-A2. Methods for reticulocyte counting (automated blood cell counters, flow cytometry, and supravital dyes). Approved Guideline - Second Edition.
2. Clinical and Laboratory Standards Institute, H26-A2. Validation, verification, and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.

QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at www.streck.com.

The brand and product names of the instruments are trademarks of their respective holders.

See www.streck.com/patents for patents that may be applicable to this product.

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