CVA Calibration Verification Assessment & CVA for CELL-DYN®

INSTRUCTIONS FOR USE INTENDED USE

CVA (Calibration Verification Assessment) and CVA for CELL-DYN® are assayed linearity control kits which can be used to assess calibration and to verify patient reportable ranges of three and five-part hematology instruments as well as automated and semi-automated non-differential hematology analyzers.

SUMMARY AND PRINCIPLES

CAP requirements¹ and CLIA regulations² both mandate that laboratories substantiate their test methods throughout the reportable range for patient test results upon installation of an analyzer. Other circumstances which may call for calibration verification include major preventive maintenance, unusual trends in control performance, or whenever recommended by the instrument manufacturer. CVA contains WBC, RBC/Hgb and Platelet concentrations which span typical patient reportable ranges; allowing the user to comply with these guidelines. CVA is vialled as separate components enabling calibration assessment without matrix interference.

Assay values are derived from replicate testing on instruments operated and maintained according to manufacturer's instructions.

REAGENTS

CVA and CVA for CELL-DYN contain stabilized human red and white blood cells and simulated or human platelets vialled as separate components in buffered aqueous media.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- CVA is not to be used for instrument calibration.
- SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

CVA/CVA for CELL-DYN is stable through the expiration date when stored at 2 °C to 10 °C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain assay values may indicate product deterioration. To determine the source of potential problems:

- Check the expiration date and discard outdated product. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Gross hemolysis (darkly colored supernatant) may indicate product deterioration. However, moderately colored supernatant is normal and should not be mistaken for product deterioration.
- Review the operating procedures of your instrument.
- If problems persist, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

Preparation for CVA and CVA for CELL-DYN

- Perform routine maintenance, quality control procedures, calibration and precision checks according to your laboratory protocol. If acceptable, proceed with running CVA or CVA for CELL-DYN.
- Remove vials of control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
- Ensure that sufficient, unexpired reagents are available and that no patient samples will be run on the instrument while running CVA or CVA for CELL-DYN.

Mixing and Handling for CVA & CVA for CELL-DYN

- To mix: **(Do NOT mix mechanically or vortex unless otherwise directed.)**
For a video demonstration, visit streck.com/mixing.
 - Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.
 - Mix vigorously by holding the vial by the ends between the thumb and finger and rapidly inverting the vial 20 times end-over-end, using a very quick turning motion of the wrist. **Examine** the bottom of the vial to make sure no cell button is visible; repeat this mixing step as needed to dislodge any cells from the bottom of the vial.
 - Continue** to invert the vial end over end 20 times, or until all cells are in suspension, and no visible cell aggregates are observed in the vial. Mix the vials in between each sample analysis by inverting in this manner 8-10 times.
 - Vials with extremely high concentrations of cells (PLT and WBC in particular, Levels 4 and beyond) may require extra mixing. Repeat steps a-c as necessary.

Important Mixing Note for CVA for CELL-DYN Platelet vials

- The **CVA for CELL-DYN Platelet vials** require an aggressive vortex mix.
For a video demonstration, visit streck.com/mixing.
 - Begin by holding the vial vertically and rolling each vial between the palms of the hands for 15-20 seconds.
 - Mix vigorously by holding the vial by the ends between the thumb and finger and rapidly inverting the vial 20 times end-over-end, using a very quick turning motion of the wrist. **Examine** the bottom of the vial to make sure no cell button is visible; repeat this mixing step as needed to dislodge any cells from the bottom of the vial.
 - Vortex the vial for a full 30 seconds.**
 - Continue** to invert the vial end over end 20 times, or until all cells are in suspension, and no visible cell aggregates are observed in the vial. Mix the vials in between each sample analysis by inverting in this manner 8-10 times.

Important Notes for CVA

- Siemens Healthcare Diagnostics ADVIA® 120/2120 Hematology System – Refer to WBCB for white blood cell count.
- Beckman Coulter® instruments – Disable the differential and blood detector (if applicable) before running CVA. (CVA may be run in primary or secondary mode.)
- Boule Medical AB MEDONIC instruments – Disable the blood detector (if applicable) before running CVA.

Important Note for CVA for CELL-DYN

- For the Emerald instrument, run CVA for CELL-DYN in the instrument's Linearity Mode.

Dilutions

- Materials necessary for dilutions but not supplied.
 - Plain tubes compatible with instrument aspiration, open or closed mode.
 - Laboratory saline or instrument-specific diluent.
 - Precision pipets, various volumes; see instrument-specific requirements.
- Dilutions - Assay: Ensure that the appropriate sample dilutions have been performed as stated on the assay sheet.
- OPTIONAL - Dilutions - All instruments: Customers desiring a linearity graph that extends lower than the assay values can perform a serial dilution according to the dilution ratios listed on the enclosed worksheet.

Analysis Protocol for CVA and CVA for CELL-DYN

- Perform an instrument background count per the instrument manufacturer's directions. If the background is within specified limits, proceed.
- RBC and Hgb results will be obtained from the RBC vials. WBC results will be obtained from the WBC vials and PLT will be obtained from the PLT vials. Instruments may report values for other parameters, but do not need to be recorded or reported.
- Refer to the assay sheet for dilution preparation or special instructions for the following instruments: Abbott CELL-DYN® 3500, 3700, 3200, Ruby and Emerald, Beckman Coulter A⁶-T⁵ 5diff, LH 500 and Mindray BC-3200.**
- Analyze all levels from one component set in order of lowest concentration to highest concentration. Run samples in Patient Mode, unless otherwise noted.
- Run each vial four times, mixing before each aspiration as indicated in mixing and handling step c. Record data on data sheet and calculate the average (mean) values.
- Repeat steps 4-5 for remaining components.
- After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both vial and cap before replacing cap and returning to refrigeration.
- Instruments with assay-specific values should verify that the mean values recovered are within the expected ranges.
 - If the mean values are outside of your laboratory's established acceptance limits, rerun CVA.
 - If the mean values remain outside of your established limits, contact your instrument service representative. Do not adjust your instrument.
- Key data into the product-specific **STATS®** Data Input Form located at streck.com and email the completed form to statsdata@streck.com, or login to your **STATS-Link®** account to submit your data.

LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with this product. The white blood cells simulate white blood cells in size, not morphology.

EXPECTED RESULTS

Values recovered should fall within the suggested tolerance limits on the accompanying assay. Due to lab-to-lab variation, some instruments may not recover values within the stated assay ranges for levels that are exceptionally concentrated or dilute. Instrument calibration, reagents and operator technique may yield differences in results obtained. Also, each laboratory must evaluate the data and use it as a guide in determining its level of instrument or testing accuracy, as well as evaluating acceptable reportable ranges. Performance goals for linearity should also be established by a laboratory based on clinical needs and what is technically achievable with instrumentation. Limits should not exceed the instrument manufacturer's linearity specification stated in the instrument manual.

- If the mean values are outside of your laboratory's established acceptance limits, rerun CVA.
- If the mean values remain outside your established limits, contact your instrument service representative. **DO NOT ADJUST YOUR INSTRUMENT.**

The following reference methods are used to whole blood calibrate instruments prior to system specific value assignment as indicated in CLSI H26-A2.³

WBC - A series of 1:500 dilutions are made using class A glassware. The lytic agent is placed in the initial dilution flask before diluting to volume. The diluting agent is an isotonic non-interfering solution. The samples are counted on a single aperture volume displacing particle counter.

RBC - A series of 1:50,000 dilutions are made using class A glassware. The diluting agent is an isotonic non-interfering solution. The samples are counted on a single aperture volume displacing particle counter.

Hb - Hemoglobin concentration is determined by converting hemoglobin to hemiglobincyanide (HiCN) and measuring absorbance at 540nm according to CLSI H15-A34 and ICSH recommendations. Hemoglobin concentration is calculated using millimolar absorption of 11.0.

Hct - Microhematocrit values are done in replicate on each sample, with capillary tubes filled and centrifuged according to the CLSI H07-A35 document. K₂EDTA is used as the anticoagulant for the collection of fresh specimens. The packed cell volume (hematocrit) is read directly using a precision metric scale. No correction is made for trapped plasma.

Plt - Platelet counting by the RBC/platelet ratio method (ICSH) is utilized as the reference method.

REFERENCES

- College of American Pathologists, 1994, Commission on Laboratory Accreditation Inspection Checklist HEM.25760.
- Clinical Laboratory Improvement Amendments of 1988, final rule: 42 CFR 493.1253.
- Clinical Laboratory Standards Institute, H26-A2, Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
- Clinical Laboratory Standards Institute, H15-A3, Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard - Third Edition.
- Clinical Laboratory Standards Institute, H7-A3, Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

The brand and product names of the instruments are trademarks of their respective holders.

See streck.com/patents for patents that may be applicable to this product.

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RÉFÉRENCES

1. College of American Pathologists, 1994, Commission on Laboratory Accreditation Inspection Checklist HEM.25760.
2. Clinical Laboratory Improvement Amendments of 1988, final rule: 42 CFR 493.1253.
3. Clinical Laboratory Standards Institute, H26-A2, Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
4. Clinical Laboratory Standards Institute, H15-A3, Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard - Third Edition.
5. Clinical Laboratory Standards Institute, H7-A3, Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

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QUELLENANGABEN

1. College of American Pathologists, 1994, Commission on Laboratory Accreditation Inspection Checklist HEM.25760.
2. Clinical Laboratory Improvement Amendments of 1988, final rule: 42 CFR 493.1253.
3. Clinical Laboratory Standards Institute, H26-A2, Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
4. Clinical Laboratory Standards Institute, H15-A3, Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard -Third Edition.
5. Clinical Laboratory Standards Institute, H7-A3, Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

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- College of American Pathologists, 1994, Commission on Laboratory Accreditation Inspection Checklist HEM.25760.
- Clinical Laboratory Improvement Amendments of 1988, final rule: 42 CFR 493.1253.
- Clinical Laboratory Standards Institute, H26-A2, Validation, verification and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.
- Clinical Laboratory Standards Institute, H15-A3, Reference and selected procedures for the quantitative determination of hemoglobin in blood. Approved Standard -Third Edition.
- Clinical Laboratory Standards Institute, H7-A3, Procedure for determining packed cell volume by the microhematocrit method. Approved Standard - Third Edition.

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Trombocyter -

REFERENSER

1. College of American Pathologists, 1994, Commission on Laboratory Accreditation Inspection Checklist HEM.25760.
2. Clinical Laboratory Improvement Amendments of 1988, final rule: 42 CFR 493.1253.
3. Clinical Laboratory Standards Institute. Validation, Verification and Quality Assurance of Automated Hematology Analyzers- Second Edition. Approved Standard.
4. Clinical Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood-Third Edition. Approved Standard.
5. Clinical Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhematocrit Method- Third Edition. Approved Standard.

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