## INSTRUCTIONS FOR USE

### INTENDED USE

MDx-Chex for BCID2 is intended for use as an external positive and negative assayed control to monitor the performance of the qualitative detection of yeast, Gram-positive and Gram-negative bacteria, as well as associated antimicrobial resistance genes, by the BIOFIRE® Blood Culture Identification 2 (BCID2) Panel on BIOFIRE FilmArray® systems. **Control 1-GN:** Gram-negative bacteria: *Acinetobacter calcoaceticus-baumannii* complex, *Bacteroides fragilis*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae* group, *Proteus* spp., *Salmonella* spp., *Serratia marcescens*, *Haemophilus influenzae*, *Neisseria meningitidis*, *Pseudomonas aeruginosa*, *Stenotrophomonas maltophilia*; antimicrobial resistance genes: KPC, CTX-M, IMP, NDM, OXA-48-like, VIM, *mcr-1*. **Control 2-GPY:** Gram-positive bacteria: *Enterococcus faecalis*, *Enterococcus faecium*, *Listeria monocytogenes*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Staphylococcus lugdunensis*, *Streptococcus agalactiae*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*; yeast: *Candida albicans*, *Candida auris*, *Candida glabrata*, *Candida krusei*, *Candida parapsilosis*, *Candida tropicalis*, *Cryptococcus neoformans/gatti*; antimicrobial resistance genes: *mecA/C* and MREJ, *vanA/B*. This product is not intended to replace manufacturer controls provided with the device.

### SUMMARY AND PRINCIPLES

Sepsis (defined as system inflammatory response syndrome in response to infection) is the 11th leading cause of death in the United States<sup>1</sup>. Life-threatening bacterial and fungal sepsis currently strikes approximately 240 out of 100,000 people per year in the U.S. (750,000 total cases), with severe sepsis (associated with acute organ dysfunction) in 95 out of 100,000 people<sup>2</sup>. Timely diagnosis and administration of effective treatment can significantly reduce mortality, duration of hospital stays, and costs due to sepsis. The BIOFIRE Blood Culture Identification 2 (BCID2) Panel simultaneously tests a single positive blood culture sample to provide results for 33 different organisms and organism groups that cause bloodstream infections and 10 genes that are known to confer antimicrobial resistance. The test can be performed on blood culture bottles that are flagged as positive by a continuously monitoring blood culture instrument. BIOFIRE BCID2 Panel results are available within approximately one hour. Rapid identification of the organism(s) in the blood culture, along with information about antimicrobial resistance gene status for select microorganisms, may aid the physician in making appropriate treatment decisions.

MDx-Chex for BCID2 is a quality control containing stabilized blood components, blood culture media components, and inactivated microorganisms resulting in a full-process, cellular-based control for the BIOFIRE BCID2 Panel. Use of full-process cellular controls are necessary to evaluate the entire analytical process, including sample lysis, nucleic acid isolation and purification, amplification, detection, and analysis, as well as the impact of PCR inhibitors and preanalytical variables. Routine use of full process quality controls can help identify variations in the test system that can lead to incorrect results.

### REAGENTS

MDx-Chex for BCID2 contains stabilized human leukocytes and erythrocytes, and the following inactivated bacteria (see Table 1) in a simulated blood culture media. Each Control tube contains enough reagent for 5 tests.

### PRECAUTIONS

- MDx-Chex for BCID2 is for *In Vitro* Diagnostic Use.
- CAUTION:** All blood products should be treated as potentially infectious. All human source material used to manufacture this product was previously established to be negative for the target analytes by a third party; non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA, and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS), West Nile Virus and Chagas disease. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
- CAUTION:** All bacterial products should be treated as potentially infectious. Source material from which this product was derived was inactivated and tested in accordance with CDC/USDA "Guidance on the Inactivation or Removal of Select Agents and Toxins for Future Use." These procedures cannot offer assurance that products containing bacteria are non-infectious.
- This product should not be disposed in general waste but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates the use of the product.

### STORAGE

MDx-Chex for BCID2 is stored at 2 °C to 25 °C. Once opened, the product may be used until the expiration date.

### INDICATION OF PRODUCT DETERIORATION

Discoloration of the product may be caused by overheating or freezing during shipping or storage. Dark colored (gross hemolysis) supernatant may be indicative of product deterioration. However, light colored (moderate hemolysis) or cloudy supernatant is normal and should not be confused with deterioration of the product.

### INSTRUCTIONS FOR USE

- If refrigerated, remove vials of control from the refrigerator and allow to sit at room temperature to acclimate for 15 minutes before use.  
Note: Always use aseptic technique when handling control tubes to prevent cross-contamination or environmental contamination.
- Prepare and hydrate a BCID2 pouch according to BIOFIRE BCID2 Panel Instructions for Use<sup>3</sup>.
- Immediately prior to use, vortex the control sample for 30 seconds to mix.  
Note: Verify the product has been adequately mixed by inverting the vial and examining the bottom for the absence of cellular material.
- Flash spin the control to remove material from cap.
- Mix the product via pipet (or transfer pipet) by repeatedly pipetting up and down 5-10 times.
- Withdraw 200µl from the control vial and add the sample to the FilmArray Sample Buffer. Close sample injection vial and invert 10 times to mix.
- Inject sample into BCID2 pouch and analyze control sample as described for patient samples in BIOFIRE's BCID2 Panel Instructions for Use<sup>3</sup>.
- Store remaining material at 2 °C to 25 °C for subsequent use.

Table 1: MDx-Chex for BCID2 Control 1-GN and Control 2-GPY Result Summary

Antimicrobial Resistance Genes		
Gene	Control 1-GN	Control 2-GPY
CTX-M	Detected	Not Detected
IMP	Detected	Not Detected
KPC	Detected	Not Detected
<i>mcr-1</i>	Detected	Not Detected
<i>mecA/C</i>	Not Detected	Detected
<i>mecA/C</i> and MREJ (MRSA)	Not Detected	Detected
NDM	Detected	Not Detected
OXA-48-like	Detected	Not Detected
<i>vanA/B</i>	Not Detected	Detected
VIM	Detected	Not Detected

Gram-Positive Bacteria		
Pathogen	Control 1-GN	Control 2-GPY
<i>Enterococcus faecalis</i>	Not Detected	Detected
<i>Enterococcus faecium</i>	Not Detected	Detected
<i>Listeria monocytogenes</i>	Not Detected	Detected
<i>Staphylococcus</i> spp.	Not Detected	Detected
<i>Staphylococcus aureus</i>	Not Detected	Detected
<i>Staphylococcus epidermidis</i>	Not Detected	Detected
<i>Staphylococcus lugdunensis</i>	Not Detected	Detected
<i>Streptococcus</i> spp.	Not Detected	Detected
<i>Streptococcus agalactiae</i> (Group B)	Not Detected	Detected
<i>Streptococcus pneumoniae</i>	Not Detected	Detected
<i>Streptococcus pyogenes</i> (Group A)	Not Detected	Detected

Gram-Negative Bacteria		
Pathogen	Control 1-GN	Control 2-GPY
<i>Acinetobacter calcoaceticus-baumannii</i> complex	Detected	Not Detected
<i>Bacteroides fragilis</i>	Detected	Not Detected
<i>Enterobacterales</i>	Detected	Not Detected
<i>Enterobacter cloacae</i> complex	Detected	Not Detected
<i>Escherichia coli</i>	Detected	Not Detected
<i>Klebsiella aerogenes</i>	Detected	Not Detected
<i>Klebsiella oxytoca</i>	Detected	Not Detected
<i>Klebsiella pneumoniae</i> group	Detected	Not Detected
<i>Proteus</i> spp.	Detected	Not Detected
<i>Salmonella</i> spp.	Detected	Not Detected
<i>Serratia marcescens</i>	Detected	Not Detected
<i>Haemophilus influenzae</i>	Detected	Not Detected
<i>Neisseria meningitidis</i>	Detected	Not Detected
<i>Pseudomonas aeruginosa</i>	Detected	Not Detected
<i>Stenotrophomonas maltophilia</i>	Detected	Not Detected

Yeasts		
Pathogen	Control 1-GN	Control 2-GPY
<i>Candida albicans</i>	Not Detected	Detected
<i>Candida auris</i>	Not Detected	Detected
<i>Candida glabrata</i>	Not Detected	Detected
<i>Candida krusei</i>	Not Detected	Detected
<i>Candida parapsilosis</i>	Not Detected	Detected
<i>Candida tropicalis</i>	Not Detected	Detected
<i>Cryptococcus neoformans/gatti</i>	Not Detected	Detected

**LIMITATIONS**

MDx-Chex for BCID2 is to be used for the BIOFIRE BCID2 Panel, on the BIOFIRE FilmArray 2.0 or Torch Systems only. It is not intended for controlling other tests or procedures.

Quality control materials should be used in accordance with local, state, federal regulations, and accreditation requirements.

**EXPECTED RESULTS**

When analyzed by the BIOFIRE BCID2 Panel, all organisms and resistance genes stated in the control should be "Detected" or "Not Detected," as indicated (see Table 1).

**REPRESENTATIVE PERFORMANCE DATA:**

Three lots of MDx-Chex for BCID2 kits were tested at 4 clinical testing sites using the BIOFIRE BCID2 Panel assay on BIOFIRE FilmArray Torch and 2.0 systems. Testing incorporated three different BCID2 pouch lots and multiple operators and instruments. A total of 240 external controls were tested over 10 days at each site. Ten Positive controls gave initial false negative results, and 1 Negative control gave an initial false positive result. All runs produced correct results upon a single retest. The overall Positive Percent Agreement (PPA) was 95.8% and Negative Percent Agreement (NPA) was 99.6%.

Summary of Multisite Precision (Reproducibility) Results									
MDx-Chex for BCID2 <sup>1</sup>	Site #1		Site #2		Site #3		Site #4		Percent Agreement (All sites combined)
	# Observed / # Expected	Percent Agreement	# Observed / # Expected	Percent Agreement	# Observed / # Expected	Percent Agreement	# Observed / # Expected	Percent Agreement	
<b>Positive:</b> Control 1-GN and Control 2-GPY, Combined	59/60	98.3%	57/60	95%	58/60	96.7%	56/60	93.3%	95.8% (230/240)
<b>Negative:</b> Control 1-GN and Control 2-GPY, Combined	60/60	100%	60/60	100%	60/60	100%	59/60	98.3%	99.6% (239/240)

<sup>1</sup> Expected result for the Positive Control is positive. Expected result for the Negative Control is negative. Denominator = total # of results for Control 1-GN and Control 2-GPY

**REFERENCES**

1. National Vital Statistics Reports, Deaths: Preliminary Data for 2010. Available from: [http://www.cdc.gov/nchs/data/nvsr/nvsr60/nvsr60\\_04.pdf](http://www.cdc.gov/nchs/data/nvsr/nvsr60/nvsr60_04.pdf).
2. Angus, D.C., et al., Epidemiology of severe sepsis in the United States: analysis of incidence, outcome, and associated costs of care. Crit Care Med, 2001. 29(7): p. 1303-10.
3. BIOFIRE Blood Culture Identification 2 (BCID2) Panel Instructions for Use - IT10048.

**ORDERING INFORMATION**

Please call our Customer Service Department 800-228-6090 for assistance. Additional information can be found online at [streck.com](http:// Streck.com).

**TECHNICAL SUPPORT**

Please call Streck Technical Services at 800-843-0912 for assistance. Additional information can be found online at [streck.com](http:// Streck.com).

**Rx Only****GLOSSARY OF SYMBOLS**

See the Instructions (IFU) tab under Resources on the product page at [streck.com](http:// Streck.com).

See [streck.com/patents](http:// Streck.com/patents) for patents that may be applicable to this product.

The brand and product names of the instruments are trademarks of their respective holders.

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## QUELLENANGABEN

1. National Vital Statistics Reports, Deaths: Preliminary Data for 2010. Verfügbar unter: [http://www.cdc.gov/nchs/data/nvsr/nvsr60/nvsr60\\_04.pdf](http://www.cdc.gov/nchs/data/nvsr/nvsr60/nvsr60_04.pdf).
2. Angus, D.C., et al., Epidemiology of severe sepsis in the United States: analysis of incidence, outcome, and associated costs of care. Crit Care Med, 2001. 29(7): S. 1303-10.
3. BIOFIRE Blood Culture Identification 2 (BCID2) Panel Instructions for Use - ITI0048.

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