INSTRUCTIONS FOR USE**INTENDED USE**

Cell-Chex® is an assayed control intended for monitoring total cell counts performed manually using a hemacytometer to validate quantitation of red and white blood cells in patient cerebrospinal fluid and body fluid samples including pleural, pericardial, peritoneal and synovial fluid. Level 1-UC contains monosodium urate crystals, and Level 1-CC contains Calcium Pyrophosphate Dihydrate (CPPD) crystals, which can be used to monitor the presence of crystals in synovial fluid.

Cell-Chex is also intended for monitoring white blood cell differentiation (Mononuclear and Polymorphonuclear; Neutrophils, Eosinophils, Basophils, Lymphocytes and Monocytes) in body fluid samples performed using Cytospin smears.

SUMMARY AND PRINCIPLES

Microscopic evaluation of cerebrospinal fluid and other body fluids is an important part of the routine analysis of body fluid specimens. Total cell count, crystal identification and white cell differential count can provide valuable information regarding infection, hemorrhage, malignancy, inflammatory conditions and irritation of joint spaces. Total cell counts and crystal identification are performed with a hemacytometer using undiluted or diluted body fluid. Crystal analysis can be performed using light or polarized light microscopy following routine lab protocol for crystal identification. Differential cell counts of body fluids are performed using Cytospin smears stained using routine hematology staining techniques.

When handled like a patient sample, Cell-Chex body fluid control will provide total cell counts and differential cell counts within the expected ranges indicated on the assay sheet. Cell-Chex is for use by laboratory professionals.

REAGENTS

Cell-Chex body fluid control contains stabilized human red blood cells, human white blood cells, monosodium urate crystals (Level 1-UC only) and CPPD crystals (Level 1-CC only) in a preservative medium. No additional materials are provided or required.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by addition of any materials to the product as supplied invalidates any diagnostic use of the product. Exception: If Level 2 control is diluted as indicated in 3.a. under Instructions for Use.
5. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

Cell-Chex is stable through the expiration date when stored unopened at 2 °C to 8 °C. After initial opening, Cell-Chex is stable for 30 days when stored at 2 °C to 8 °C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values for total cell counts or the inability to differentiate the cell types may indicate product deterioration. If the recovered values are not within the expected ranges:

1. Review the control product package insert Instructions For Use and Limitations sections.
2. Check the expiration date of Cell-Chex. Discard outdated product.
3. Assay an unopened vial of Cell-Chex. If the values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

1. Remove the controls from the refrigerator. It is not necessary to warm the controls to room temperature before using.
2. To mix: **(Do NOT mix mechanically or vortex.)**
For a video demonstration, visit streck.com/mixing.
 - a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.
 - d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period.
3. Samples must be removed using a clean capillary tube or pipette tip. The vial must be closed immediately after sampling is complete. Care should be taken to prevent cross-contamination of the control.
 - a. For hemacytometer counts, refer to your hemacytometer package insert for specific instructions. The Level 2 control is concentrated enough to perform a dilution and obtain a reportable value.
 - b. If using Cell-Chex for synovial fluid control and crystal analysis, further crystal evaluation of a positive control can be performed using polarized light microscopy or by following established lab protocol for synovial fluid analysis.
 - c. For white blood cell differential counts, prepare a Cytospin smear according to your established

laboratory procedure. For best results, use 250µl of control product. **Addition of any material to Cell-Chex prior to preparation of a smear is not necessary or recommended.** Stain the Cytospin smear of Cell-Chex according to your established laboratory procedure for routine hematology. **The differential preparation used to obtain Cell-Chex assay values is centrifuged for 7 minutes at 500 RPMs.**

4. After sampling, return to refrigeration for maximum open-vial stability. Wipe the threads of both the vial and cap before replacing cap and returning to refrigeration. Report any serious incident to Streck and appropriate regulatory entities including the competent authority of the member state in which the user and/or patient is established as applicable.

LIMITATIONS

1. Analysis of Cell-Chex cannot be performed on an automated hematology analyzer. The stabilization process of the red blood cells will cause them to interfere with the white blood cell count.
2. Cell-Chex white blood cells are stabilized. Pretreatment of Cell-Chex prior to preparation of a Cytospin smear is not necessary.

EXPECTED RESULTS AND THEIR DERIVATION

The assigned values for red blood cells, white blood cells and differential values were determined by replicate analyses of hemacytometer counts and stained Cytospin smears. The assigned values of the total cell counts are based on the average number of cells per µL counted in the nine large gridded squares on both sides of the Neubauer Improved hemacytometer counting chamber. The assigned values of the differential cell counts are based on replicate 100-cell differentials. Assay range values were established based on +/- 3SD for WBC differential parameters and +/- 2SD for the RBC and WBC parameters. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control.

PERFORMANCE CHARACTERISTICS

The assigned cellular component values are presented as a mean and a recovery range. The mean cellular component assay values are derived from repetitive testing on multiple vials of product using Neubauer Improved hemacytometers, which are operated and maintained according to the manufacturer's instructions. White cell differential cellular assay values were determined by evaluating repetitive Cytospin smears from different vials run according to the manufacturer's instructions. Crystals are identified as "present" or "not present", do not have a recovery range, and are determined visually using hemacytometers on representative vials. The recovery ranges for cellular components are intended to reflect laboratory variability, thus, they are within +/- 3SD for WBC differential parameters and +/- 2SD for the RBC and WBC parameters. Always perform quality control according to good laboratory practice, laboratory director's requirements, and any regulatory or accredited requirements.

QUALITY CONTROL PROGRAM

Streck offers STATS®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only**GLOSSARY OF SYMBOLS**

See the Instructions (IFU) tab under Resources on the product page at streck.com.

CHANGES FROM PREVIOUS VERSION

Updated to meet requirements of Regulation (EU) 2017/746 and EU Authorized Representative address.

See streck.com/patents for patents that may be applicable to this product.

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