INSTRUCTIONS FOR USE**INTENDED USE**

Cyto-Chex® BCT is intended for collection and storage of blood specimens for immunophenotyping of white blood cells by flow cytometry. Recovery of lymphocyte subset cell markers of the HIV panel can be accomplished over a 14-day period following collection.

SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes in a blood sample. Logistics, such as transport of samples from collection site to analysis location, create situations in which sample analysis cannot be performed immediately after collection. Sample preservation for future analysis provides a solution which addresses such situations.

Cyto-Chex BCT is a direct-draw blood collection tube containing an anticoagulant and preservative, which preserves the cell surface antigens of white blood cells (leukocytes) until processing and analysis can be performed.

Subsets of leukocytes can be distinguished on the basis of cell surface antigens using fluorescent antibodies and flow cytometry.^{1,2} Qualitative and quantitative changes in leukocyte subsets are used to identify and monitor immunodeficiency and hematologic diseases.

Cyto-Chex BCT is designed to preserve peripheral blood samples' qualitative and quantitative leukocyte subset characteristics. Cyto-Chex BCT is for use by laboratory professionals.

REAGENTS

Cyto-Chex BCT contains the anticoagulant, K₃EDTA, and a cell preservative in a liquid medium. No additional materials are provided or required.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- Do not freeze specimens collected in Cyto-Chex BCT as breakage could result.
- Do not use tubes after expiration date.
- Do not use tubes for collection of materials to be injected into patients.
- Product is intended for use as supplied. Do not dilute or add other components to Cyto-Chex BCT.
- Overfilling or under-filling of tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.

CAUTION

- Glass has the potential for breakage; precautionary measures should be taken during handling.
 - All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection. Dispose of in accordance with federal, state and local regulations. Avoid contact with skin and mucous membranes.
 - Product should be disposed with infectious medical waste.
 - Remove stopper by either gently rocking the stopper from side to side or by grasping with a simultaneous twisting and pulling action. A "thumb roll" procedure for stopper removal is NOT recommended, as tube breakage and injury may result. Reinsert stopper by gently pushing stopper onto tube with a simultaneous twisting action.
 - A pediatric tube adapter is strongly recommended to better fit the 10.25mm diameter **2.0ml** Cyto-Chex BCT.
- SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

- When stored at 2 °C to 30 °C, empty Cyto-Chex BCT is stable through expiration date.
- Short term storage at 2 °C to 40 °C is acceptable for empty Cyto-Chex BCT for up to 14 days.
- Do not freeze Cyto-Chex BCT. Proper insulation may be required for shipment during extreme temperature conditions.
- HIV Panel of Markers**
 - The HIV panel of markers including CD3, CD4, CD8, CD16/CD56, CD19 and CD45 are stable in blood samples stored in Cyto-Chex BCT for up to 14 days at 18 °C to 22 °C and up to 7 days at 30 °C.
 - CD8 is stable for up to 48 hours at 37 °C.
 - CD3 and CD4 are stable for up to 72 hours at 37 °C.


INDICATIONS OF PRODUCT DETERIORATION

- Cloudiness or precipitate visible in reagent of empty tube.
- If indications of product deterioration occur, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

- Collect specimen by venipuncture according to CLSI PRE02.³

Note: For the 2.0ml Cyto-Chex BCT, care must be taken to center the tube so the middle of the stopper is punctured per the figure below. To aid with proper insertion of the tube in standard holders and stabilize the tube during the draw, a pediatric tube adapter is strongly recommended to modify the standard holder to fit the 10.25mm diameter **2.0ml** Cyto-Chex BCT (Streck part number: 213559, 218980).


Correctly punctured stopper

Prevention of Backflow – Since Cyto-Chex BCT contains chemical additives, it is important to avoid possible backflow from the tube.

To guard against backflow, observe the following precautions:

- Keep patient's arm in the downward position during the collection procedure.
 - Hold the tube with the stopper in the uppermost position so that the tube contents do not touch the stopper or the end of the needle during sample collection.
 - Release tourniquet once blood starts to flow in the tube, or within 2 minutes of application.
- Follow recommendations for order of draw outlined in CLSI PRE02.³ Cyto-Chex BCT can be drawn after the EDTA tube.
 - Fill tube completely.

- Remove tube from adapter and immediately mix by gentle inversion 8 to 10 times. Inadequate or delayed mixing may result in inaccurate test results. One inversion is a complete turn of the wrist, 180 degrees, and back per the figure below:

- After collection, transport and store tubes within the recommended temperature range.
- Perform immunophenotyping by flow cytometry in accordance with instrument manufacturer's instructions.

Report any serious incident to Streck and appropriate regulatory entities including the competent authority of the member state in which the user and/or patient is established as applicable.

Note:

- Light scatter positions of cells stabilized with Cyto-Chex BCT may differ slightly from those of untreated cells.
- Cyto-Chex BCT does not dilute blood samples; therefore, no dilution factor correction is necessary to obtain absolute count values.
- As in the case with most clinical laboratory specimens, hemolysis, icteris and lipemia may affect the results obtained on blood samples preserved with Cyto-Chex BCT. The specimen condition should be noted on the laboratory report since the results may be suspect. Grossly hemolyzed samples should be rejected.

EXPECTED RESULTS AND THEIR DERIVATION

Immunophenotypic analysis performed on peripheral blood samples stored in Cyto-Chex BCT will provide the same results as when performed on fresh specimens.

CD4 + Cell Count Recovery

	EDTA tube	Cyto-Chex BCT	
Healthy Donors	Initial Cells/μl	Day 7 Cells/μl	Day 14 Cells/μl
1	597	684	641
2	1275	1421	1267
3	559	567	563
4	646	591	572
5	843	844	915
HIV+ Patients	Initial Cells/μl	Day 7 Cells/μl	Day 14 Cells/μl
1	709	760	756
2	511	561	565
3	285	299	300
4	327	340	373
5	83	73	89

LIMITATIONS

- For single use only.
- DO NOT** transfer samples into Cyto-Chex BCT that have been collected in another anticoagulant or preservative.

REFERENCES

- CDC guidelines for performing single-platform absolute CD4+ T-cell determinations with CD45 gating for persons infected with human immunodeficiency virus. MMWR 52(RR-2) 2003.
- Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
- Clinical and Laboratory Standards Institute, PRE02, Collection of diagnostic venous blood specimens. Approved Standard - Eighth Edition.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only**GLOSSARY OF SYMBOLS**

See the Instructions (IFU) tab under Resources on the product page at streck.com.

CHANGES FROM PREVIOUS VERSION

Updated to add Czech translation.

See streck.com/patents for patents that may be applicable to this product.

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LITTERATUR

1. CDC guidelines for performing single-platform absolute CD4+ T-cell determinations with CD45 gating for persons infected with human immunodeficiency virus. MMWR 52(RR-2) 2003.
2. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
3. Clinical and Laboratory Standards Institute, PRE02, Collection of diagnostic venous blood specimens. Approved Standard - Eighth Edition.

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- TÀI LIỆU THAM KHẢO**
- CDC guidelines for performing single-platform absolute CD4+ T-cell determinations with CD45 gating for persons infected with human immunodeficiency virus. MMWR 52(RR-2) 2003.
 - Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
 - Clinical and Laboratory Standards Institute, PRE02, Collection of diagnostic venous blood specimens. Approved Standard - Eighth Edition.

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