INSTRUCTIONS FOR USE**INTENDED USE**

CD-Chex Plus® is intended to be used as a quality control material for evaluating monoclonal antibody binding by flow cytometry. CD-Chex Plus control cells possess surface antigens detectable with monoclonal antibodies. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry they provide a reference for normal peripheral blood leukocytes. CD-Chex Plus is designed for use on the BD® Biosciences and Beckman Coulter® flow cytometry systems.

SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes present in a whole blood sample. Quantitation of peripheral blood leukocyte subsets has been shown to be useful in the evaluation of immunodeficiency, autoimmune, and oncologic disorders. CD-Chex Plus is a whole blood positive procedure control with assayed values used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance. CD-Chex Plus control cells mimic whole blood samples by possessing surface antigens detectable with fluorescent monoclonal antibodies by flow cytometry. CD-Chex Plus is for use by laboratory professionals.

Subsets of leukocytes can be distinguished on the basis of light scatter properties as well as the presence of cell surface antigens. CD-Chex Plus is designed to represent normal peripheral blood leukocytes where lymphocytes, monocytes, and granulocytes are distinguishable on the basis of light scatter properties.

REAGENTS

CD-Chex Plus contains stabilized human blood in a preservative medium. No additional materials are provided or required.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. **CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates the intended use of the product.
5. CD-Chex Plus should not be used as a calibrator.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.
7. This product contains preservatives which may result in sensitization or an allergic reaction; all users are to wear standard laboratory PPE.

STORAGE AND STABILITY

CD-Chex Plus is stable through the expiration date when stored unopened at 2 °C to 8 °C. After initial opening, CD-Chex Plus is stable for 30 days when stored at 2 °C to 8 °C.

DO NOT FREEZE.

INDICATIONS OF PRODUCT DETERIORATION

If CD-Chex Plus values are not within the expected range on the assay sheet or granulocytes show a loss of forward scatter (FSC):

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of CD-Chex Plus. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.
4. Clumping of the cell suspension indicates instability or deterioration, in which case the reagent should not be used.

INSTRUCTIONS FOR USE

1. CD-Chex Plus should only be used with the instruments described above in the INTENDED USE.
2. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
3. Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
4. **Mixing Procedure (mechanical mixing by vortex or rotator is not recommended):**

For a video demonstration, visit streck.com/mixing.

- a. Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- b. Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.

- c. Aliquot immediately after mixing.
- d. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.

Note: Vials stored for an extended period of time may require extra mixing.

5. Return control reagent to refrigeration immediately after sampling to ensure maximum open-vial stability.
6. Add recommended monoclonal antibodies according to manufacturer's instruction to each tube and mix gently.
7. Incubate according to antibody manufacturer's instructions.
8. Add recommended amount of RBC lysing agent and follow manufacturer's instructions.

9. Analyze by flow cytometry using your laboratory's established protocol. Report any serious incident to Streck and appropriate regulatory entities including the competent authority of the member state in which the user and/or patient is established as applicable.
10. For the Accellix™ platform, refer to Accellix Assay Instructions for Use and TBNK-16 (L) Autoclassification Data Sheet.

LIMITATIONS

1. Stabilized material is considered nonviable and not compatible with viability dyes and kits.

PROCEDURES

1. Instrument Procedure. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. Monoclonal Antibody Procedure. Use monoclonal antibodies according to manufacturer's instructions for patient samples.
3. RBC Lysing Procedure. Use RBC lysing agent according to manufacturer's instructions for patient samples.

EXPECTED RESULTS AND THEIR DERIVATION

The mean assay values provided for each parameter are derived from replicate analyses on properly compensated flow cytometers.^{1,2} The assay values are obtained using common flow cytometry reagents. See the IFU and assay for limitations or specific instructions for reagents.

Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different reagents, laboratory protocols, instrument calibration, maintenance, and operator technique.

REFERENCES

1. Clinical and Laboratory Standards Institute H43-A2 Clinical Flow Cytometric Analysis of Neoplastic Hematolymphoid Cells; Approved Guideline - Second Edition.
2. Clinical and Laboratory Standards Institute H42-A2 Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline - Second Edition.

QUALITY CONTROL PROGRAM

Streck offers **STATS®**, an interlaboratory quality control program, to all customers at no charge. For more information, contact the **STATS** Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only**GLOSSARY OF SYMBOLS**

See the Instructions (IFU) tab under Resources on the product page at streck.com.

CHANGES FROM PREVIOUS VERSION

Updated Instructions for Use Section.

EU SUMMARY OF SAFETY AND PERFORMANCE (SSP)

SSP is available on the EUDAMED website (ec.europa.eu/tools/eudamed/).

All product names, logos, brands, and marks are property of their respective owners.

See streck.com/patents for patents that may be applicable to this product.

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