INSTRUCTIONS FOR USE

INTENDED USE

ESR-Chex® is an assayed bi-level control for evaluating the accuracy and precision of automated and manual Erythrocyte Sedimentation Rate (ESR) methods.

SUMMARY

Clinical laboratories are required to establish quality control programs for automated and semi-automated procedures. Daily use of ESR-Chex provides assurance that the erythrocyte sedimentation rate procedure is being performed properly.

PRINCIPLES

ESR-Chex is designed to be handled in the same manner as a patient specimen. There are several factors concerning the environment and techniques that may affect the outcome of erythrocyte sedimentation rate methods. When these are properly controlled, the values recovered with ESR-Chex will be within the Expected Ranges printed on the assay. Accurate analysis of patient specimens is assured.

REAGENTS

ESR-Chex contains stabilized human red blood cells in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.
5. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

ESR-Chex is stable through the expiration date when stored at 2 °C-10 °C. After opening, ESR-Chex is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C-30 °C.

INDICATIONS OF PRODUCT DETERIORATION

Discoloration of the product may be caused by overheating or freezing during shipping or storage. Gross hemolysis (darkly colored supernatant) may be indicative of product deterioration. However, moderately colored supernatant is normal and should not be confused with deterioration of the product.

INSTRUCTIONS FOR USE


1. INSTRUMENT PROCEDURE. Follow the instrument manufacturer's instructions for ESR testing for automated methods.
2. MANUAL PROCEDURE. For manual testing, use ESR-Chex as you would a patient sample.
3. CAUTION. **DO NOT remove sodium citrate or sodium chloride from tubes before using this control.**
4. It is good laboratory practice to prepare a fresh sample of ESR-Chex daily to control the ESR procedure.

MIXING AND HANDLING INSTRUCTIONS

For a video demonstration, visit Streck.com/mixing.

Incomplete mixing can invalidate both the sample drawn and the remaining product in the vial.

1. Remove vials from refrigerator and allow them to equilibrate to room temperature (20-30 minutes).
2. To Mix:
 - a. **Vortex** unopened vials for up to 60 seconds to resuspend product prior to first time use.
 - b. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- c. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- d. Examine the bottom of the vial to ensure that the product is completely resuspended and repeat steps 2b-2c, if necessary. Allow bubbles to disperse and remix by inversion 8-10 times immediately prior to sampling.
 - e. Aliquot the sample immediately after mixing. If mixed vials sit for more than 1 minute before drawing the sample, remix the vial by repeating steps 2b-2d. Subsequent analyses during this test period may be performed by repeating steps 2b-2d.
3. Previously opened vials stored at room temperature (18 °C-30 °C) should not require additional vortexing. Mix previously opened vials by following steps 2b-2e and repeat as necessary to ensure the cells are completely resuspended.
 4. Follow the manufacturer's directions for the method in use to set up and fill the sedimentation rate tube.
 5. Wipe threads of vial and cap with clean tissue before closing. Recap the vial tightly.
 6. It is recommended that opened vials be stored at room temperature (18 °C-30 °C) for the duration of the open-vial stability to expedite resuspending the cells for subsequent analysis.

LIMITATIONS

ESR-Chex is to be used for Erythrocyte Sedimentation Rate testing only. It is not intended for controlling other hematology procedures.

EXPECTED RESULTS

The assay values are derived from replicate analysis on both automated and manual methods. Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits. However, the control means established by the laboratory should fall within the Expected Range specified for the control.

In order to achieve values within the published assay ranges, laboratories should prepare a fresh sample of ESR-Chex daily. Failure to adhere to this guideline may result in values outside the published expected ranges. It is advised that laboratories establish their own means and ranges when utilizing control material outside of these handling instructions.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2¹ it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

If established means do not fall within the Expected Range, please evaluate all mechanical and physical factors that could affect the outcome, such as handling, temperature, vibration, tube position and product expiration. Inadequate mixing of the product will affect the outcome of the test, see Handling Instructions. Sample must be taken within 60 seconds of thoroughly mixing the vial. Assay an unopened vial of ESR-Chex. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

REFERENCES

1. Clinical and Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.

QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

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