




Superior sed-rate accuracy in only 30 minutes.



ESR-Auto Plus® & ESR-100

Streck ESR-Auto Plus and ESR-100 automated sed-rate analyzers exhibit >98% correlation to the Modified Westergren method. Graphs shown display data collected in a comparison of 150 fresh patient samples on the ESR-Auto Plus versus the leading competitors' instruments.



ESR-Auto Plus



Streck's Instrument Features

- 98% correlation to Modified Westergren
- Random access eliminates sample batching
- Internal printer automatically reports results
- Barcode scanner for sample identification
- Results available in 30 minutes
- ESR-Chex sed-rate controls available
- Direct draw ESR-Vacuum Tubes preserve the integrity of the patient sample for up to 72 hours when transported and stored at 2-10°C.

High Volume Laboratories

The ESR-100 is a high capacity 100-position analyzer designed for laboratories that run over 100 samples daily. The ESR-100 offers the same accuracy as the ESR-Auto Plus and results are available in 30 minutes.



ESR-100

Streck
Clinical Laboratory Products
You Can Count On.



800.843.0912
www.streck.com