

# 10 Reasons to Choose the ESR-Auto Plus<sup>®</sup>

An automated sedimentation rate analyzer.



Over 1500  
units placed

- **>98% Correlation to Modified Westergren Method**  
Provides accurate results in only 30 minutes (15-minute prediction mode). Ensures confidence and saves time and money.
- **Internal Printer**  
Automatically prints patient results at the conclusion of the test. Eliminates repeat testing and saves valuable benchtop space.
- **Full Random Access**  
Eliminates sample batching and provides flexibility. Ten-sample capacity allows patient testing at your convenience.
- **Positive Sample Identification**  
Barcode scanner aids in patient sample identification and prevents transcription errors. Stores patient results in log file.
- **QC Data Storage**  
Stores lot-specific QC information for up to 100 runs of QC data per level.
- **LIS Compatible**  
An LIS interface prevents transcription errors and aids in patient sample identification.
- **72-Hour Patient Sample Stability**  
ESR-Vacuum Tubes preserve the integrity of the patient sample for up to 72 hours at 2-10° C. Ensures specimen stability during transit from off-site blood collection facility.
- **Closed Blood Collection System**  
Reduces exposure to potentially hazardous material. Allows direct vacuum blood collection or transfer from an EDTA (purple-top) tube.
- **Complete System**  
The ESR-10 Manual Rack is a 30-minute manual system that uses Streck ESR-Vacuum Tubes. Safety-coated tubes and ESR-Chex control are available, as well as an ESR placement program.
- **Excellent Service**  
Streck provides superior sales, customer and technical service to help you through the ordering and evaluation process. The ESR-Auto Plus has a proven track record and substantial peer group.

800.843.0912  
www.streck.com

**Streck**

Clinical Laboratory Products You Can Count On.