



Temp-Chex

Temp-Chex with Red Spirit thermometers monitor temperatures between -90 °C and 50 °C and come calibrated with a NIST-traceable certificate to ensure accuracy and compliance with regulatory requirements. They are mercury free and enclosed in a plastic case to protect the thermometer and contain the contents if broken, making it environmentally safe. Temp-Chex with Red Spirit is scaled in 0.5° increments (ultra-low thermometer features 1° gradation), allowing you to meet CAP and CLIA standards.

Features

- + Individually serial-numbered and NIST-traceable
- + Measures temperatures from -90 °C to 50 °C
- + Yellow background for easier reading
- + Enclosed in a durable plastic case to contain contents in case of breakage
- + Can be mounted magnetically or by adhesive
- + Downloadable wall chart for daily temperature monitoring available

Benefits

- + A wide temperature range allows Temp-Chex thermometers to be used to monitor most laboratory applications
- + Mercury free and non-hazardous for increased laboratory safety
- + Thermometer probe is housed in a vial of propylene glycol, which creates a temperature buffer that prevents fluctuations

Ordering information

Red Spirit	
Level	Description
-5 °C to 15 °C Refrigerator	240043
-30 °C to 1 °C Freezer	240046
15 °C to 50 °C Incubator	240049
-90 °C to 25 °C Ultra-Low Freezer	240047