

Whole Blood Assayed Hematology Control<sup>1</sup>

Rx Only

Open-vial stability 95 days<sup>2</sup>

 2026-03-31

350206-40  
2023-09

| Automated Methods <sup>3</sup>  | [LOT] L1 - 50901380 |         | [LOT] L2 - 50901381 |         |
|---------------------------------|---------------------|---------|---------------------|---------|
|                                 | $\bar{x}$           | +/-     | $\bar{x}$           | +/-     |
| Instrument <sup>4</sup>         | mm / hr             | mm / hr | mm / hr             | mm / hr |
| Vital Diagnostics Excyte 10 / M | 9                   | 8       | 86                  | 18      |

| Manual Methods <sup>5</sup>                     | [LOT] L1 - 50901380 |         | [LOT] L2 - 50901381 |         |
|-------------------------------------------------|---------------------|---------|---------------------|---------|
|                                                 | $\bar{x}$           | +/-     | $\bar{x}$           | +/-     |
| Instrument <sup>4</sup>                         | mm / hr             | mm / hr | mm / hr             | mm / hr |
| Classical Westergren <sup>a</sup>               | 5                   | 4       | 19                  | 7       |
| Modified Westergren - Globe Sedi-Rate           | 5                   | 4       | 42                  | 20      |
| Modified Westergren - Polymedco Sediplast       | 5                   | 4       | 49                  | 20      |
| Modified Westergren - Streck ESR-10 Manual Rack | 9                   | 8       | 65                  | 18      |
| Wintrobe                                        | 5                   | 4       | 20                  | 7       |

a Undiluted Westergren Method

Valeurs principales assignées / Unverdünnte Westergren-Methode / Metodo Westergren non diluito / Método Westergren sin diluir / Outspädd Westergrenmetod

- 1 Whole Blood Assayed Hematology Control
- Contrôle d'hématologie dosé « sang total » / Mit Vollblut getestete Hämatologiekontrolle / Controllo ematologico per analisi di sangue intero / Control de hematología de sangre entera analizada / Helblodsanalyserad hematologikontroll
- 2 Open-vial stability 95 days
- Stabilité en flacon ouvert 95 jours / Stabilität geöffneter Flaschen 95 tage / Stabilità della fiala aperta 95 giorni / Estabilidad de la cápsula abierta 95 días / Hållbarhet för öppen flaska 95 dagar
- 3 Automated Methods
- Méthodes automatisées / Automatische Verfahren / Metodi automatizzati / Métodos automatizados / Automatiska metoder
- 4 Instrument
- Instrument / Gerät / Strumento / Instrumento / Instrument
- 5 Manual Methods
- Méthodes manuelles / Manuelle Verfahren / Metodi manuali / Métodos manuales / Manuella metoder
- $\bar{x}$  Mean
- Moyenne / Mittelwert / Media / Media / Medelvärde
- +/- Expected Range
- ± Intervalle escompté / ± Erwartungsbereich / ± Range previsto / ± Intervalo previsto / ± Förväntat intervall