



Whole Blood Assayed Hematology Control<sup>1</sup>

Rx Only

Open-vial stability 95 days<sup>2</sup>



2027-04-27

350206-41  
2025-06



streck.com/IFU

**Automated Methods<sup>3</sup>**

	LOT L1 - 61171380		LOT L2 - 61171381	
	$\bar{x}$	+/-	$\bar{x}$	+/-
Instrument <sup>4</sup>	mm / hr	mm / hr	mm / hr	mm / hr
Vital Diagnostics Excyte 10 / M	9	8	71	18

**Manual Methods<sup>5</sup>**

	LOT L1 - 61171380		LOT L2 - 61171381	
	$\bar{x}$	+/-	$\bar{x}$	+/-
Instrument <sup>4</sup>	mm / hr	mm / hr	mm / hr	mm / hr
Classical Westergren <sup>a</sup>	5	4	17	7
Modified Westergren - Globe Sedi-Rate	5	4	37	20
Modified Westergren - Polymedco Sediplast	5	4	40	20
Wintrobe	5	4	18	7

**a Undiluted Westergren Method**

Valeurs principales assignées / Unverdünnte Westergren-Methode / Metodo Westergren non diluito / Método Westergren sin diluir / Outspädd Westergrenmetod

**1 Whole Blood Assayed Hematology Control**

Contrôle d'hématologie dosé « sang total » / Mit Vollblut getestete Hämatologiekontrolle / Controllo ematologico per analisi di sangue intero / Control de hematologia de sangre entera analizada / Helblodsanalyserad hematologikontroll

**2 Open-vial stability 95 days**

Stabilité en flacon ouvert 95 jours / Stabilität geöffneter Flaschen 95 tage / Stabilità della fiala aperta 95 giorni / Estabilidad de la cápsula abierta 95 días / Hållbarhet för öppen flaska 95 dagar

**3 Automated Methods**

Méthodes automatisées / Automatische Verfahren / Metodi automatizzati / Métodos automatizados / Automatiska metoder

**4 Instrument**

Instrument / Gerät / Strumento / Instrumento / Instrument

**5 Manual Methods**

Méthodes manuelles / Manuelle Verfahren / Metodi manuali / Métodos manuales / Manuella metoder

**$\bar{x}$  Mean**

Moyenne / Mittelwert / Media / Media / Medelvärde

**+/- Expected Range**

± Intervalle escompté / ± Erwartungsbereich / ± Range previsto / ± Intervalo previsto / ± Förväntat intervall



Streck  
7002 S. 109 Street, La Vista, NE 68128 USA



MediMark Europe Sarl  
11 rue Emile Zola  
38100 Grenoble, France