

Whole Blood Assayed Hematology Control¹

Rx Only

Open-vial stability 95 days²



2025-04-29

350206-40
2023-09

Automated Methods³

| | [LOT] L1 - 41201380 | | [LOT] L2 - 41201381 | |
|---------------------------------|---------------------|---------|---------------------|---------|
| | \bar{x} | +/- | \bar{x} | +/- |
| Instrument ⁴ | mm / hr | mm / hr | mm / hr | mm / hr |
| Vital Diagnostics Excyte 10 / M | 9 | 8 | 90 | 18 |

Manual Methods⁵

| | [LOT] L1 - 41201380 | | [LOT] L2 - 41201381 | |
|---|---------------------|---------|---------------------|---------|
| | \bar{x} | +/- | \bar{x} | +/- |
| Instrument ⁴ | mm / hr | mm / hr | mm / hr | mm / hr |
| Classical Westergren ^a | 5 | 4 | 21 | 7 |
| Modified Westergren - Globe Sedi-Rate | 5 | 4 | 47 | 20 |
| Modified Westergren - Polymedco Sediplast | 5 | 4 | 51 | 20 |
| Modified Westergren - Streck ESR-10 Manual Rack | 9 | 8 | 69 | 18 |
| Wintrobe | 5 | 4 | 23 | 7 |

a Undiluted Westergren Method

Valeurs principales assignées / Unverdünnte Westergren-Methode / Metodo Westergren non diluito / Método Westergren sin diluir / Outspädd Westergrenmetod

1 Whole Blood Assayed Hematology Control

Contrôle d'hématologie dosé « sang total » / Mit Vollblut getestete Hämatologiekontrolle / Controllo ematologico per analisi di sangue intero / Control de hematologia de sangre entera analizada / Helblodsanalyserad hematologikontroll

2 Open-vial stability 95 days

Stabilité en flacon ouvert 95 jours / Stabilität geöffneter Flaschen 95 tage / Stabilità della fiala aperta 95 giorni / Estabilidad de la cápsula abierta 95 días / Hållbarhet för öppen flaska 95 dagar

3 Automated Methods

Méthodes automatisées / Automatische Verfahren / Metodi automatizzati / Métodos automatizados / Automatiska metoder

4 Instrument

Instrument / Gerät / Strumento / Instrumento / Instrument

5 Manual Methods

Méthodes manuelles / Manuelle Verfahren / Metodi manuali / Métodos manuales / Manuella metoder

\bar{x} Mean

Moyenne / Mittelwert / Media / Media / Medelvärde

+/- Expected Range

± Intervalle escompté / ± Erwartungsbereich / ± Range previsto / ± Intervalo previsto / ± Förväntat intervall