

CD-Chex™ Daily

One vial. 51 markers. Your panel covered.



Control Marker Guide/Coverage

Phenotype	CD-Chex™ Daily*	CD-Chex Plus®	CD-Chex CD34®	CD-Chex Select®*	CD-Chex CD103 Plus®*	CD-Chex CD117® Plus*	CD-Chex TdT Plus®
CD1a							•
CD2	•	•					•
CD3	•	•					
cCD3							•
CD4	•	•					•
CD5	•	•					•
CD7	•	•					•
CD8	•	•					•
CD10	•	•					
CD11b	•	•					
CD11c	•	•					
CD13	•	•					
CD14	•	•					
CD15	•	•					
CD16	•	•					
CD16-CD56	•	•					
CD19	•	•					
CD20	•	•					
CD22	•	•					
CD23	•	•					
CD24	•						
CD25						•	
CD26	•			•			
CD27	•						
CD28	•						
CD30					•		
CD33	•	•					
CD34 (Progenitor)	•	•	•				
CD34 (Abnormal)							•
CD36	•						
CD38+ / CD45+	•	•				•	
CD38+ / CD45-						•	
CD41	•			•			
CD43	•						
CD45RA	•						
CD45RO	•						
CD45	•	•					
CD56	•	•			•		
CD57	•			•			
CD58	•			•			
CD61	•			•			
CD64	•			•			
CD71						•	
CD79a	•			•			
CD79b	•						
CD81	•						
CD103					•		
CD117						•	
CD138					•		
CD200	•						
FMC7	•			•			
HLA-DR	•	•					
IgM	•						
Kappa / CD19	•	•					
Lambda / CD19	•	•					
cKappa	•						
cLambda	•				•		
cMPO	•			•			
TCR α/β	•			•			
TCR γ/δ	•			•			
TdT							•
TRBC1	•						
Lymphocyte	•	•					
Monocyte	•	•					
Granulocyte	•	•					
WBC	•	•	•	•	•	•	•

• CD-Chex Plus RUO markers will be discontinued with lot 6348, shipping Dec. 14, 2026.

• CD-Chex Select will be discontinued with lot 6208, shipping July 27, 2026.

*CD-Chex Daily, CD-Chex Select, CD-Chex CD103 Plus, CD-Chex CD117 Plus, and CD-Chex TdT Plus are For Research Use Only. Not for use in diagnostic procedures.