




Quality Control for Multiplex Molecular Respiratory Assays

350868  
2026-06

For Research Use Only.  
Not for use in diagnostic procedures.

 2027-06-07

The product may be used until expiration date.

**Expected Results**

Viruses			
Target	LOT 61591230	Positive Control	LOT 61591231 Negative Control
Coronavirus (Seasonal): HKU1, NL63, OC43, 229E		Detected	Not Detected
Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2)		Detected	Not Detected
Human Metapneumovirus		Detected	Not Detected
Human Rhinovirus/ Enterovirus		Detected	Not Detected
Influenza A		Detected	Not Detected
Influenza A H1-2009		Detected	Not Detected
Influenza A H3		Detected	Not Detected
Influenza B		Detected	Not Detected
Parainfluenza Virus 1		Detected	Not Detected
Parainfluenza Virus 2		Detected	Not Detected
Parainfluenza Virus 3		Detected	Not Detected
Parainfluenza Virus 4		Detected	Not Detected
Respiratory Syncytial Virus		Detected	Not Detected
Bacteria			
<i>Bordetella parapertussis</i> (IS1001)		Detected	Not Detected
<i>Bordetella pertussis</i> (ptxP)		Detected	Not Detected
<i>Chlamydia pneumoniae</i>		Detected	Not Detected
<i>Mycoplasma pneumoniae</i>		Detected	Not Detected
<i>Streptococcus dysgalactiae</i> (Group C/G Strep)		Detected	Not Detected
<i>Streptococcus pyogenes</i> (Group A Strep)		Detected	Not Detected

All product names, logos, brands, and marks are property of their respective holders.