

Quality Control for the Luminex VERIGENE[®] Gram-Negative Blood Culture Test (BC-GN)¹

350786-2
2025-10

2027-09-21

Rx Only

The product may be used until expiration date².

Luminex VERIGENE[®] System

Gram-Negative Bacteria ³			
Target ⁴	LOT 62641159	Positive Control ⁵	LOT 62641156 Negative Control ⁶
<i>Escherichia coli</i>		Detected ⁷	Not Detected ⁸
<i>Klebsiella pneumoniae</i>		Detected ⁷	Not Detected ⁸
<i>Klebsiella oxytoca</i>		Detected ⁷	Not Detected ⁸
<i>Pseudomonas aeruginosa</i>		Detected ⁷	Not Detected ⁸
<i>Acinetobacter</i> spp.		Detected ⁷	Not Detected ⁸
<i>Citrobacter</i> spp.		Detected ⁷	Not Detected ⁸
<i>Enterobacter</i> spp.		Detected ⁷	Not Detected ⁸
<i>Proteus</i> spp.		Detected ⁷	Not Detected ⁸

Antimicrobial Resistance Genes ⁹			
Gene ¹⁰	LOT 62641159	Positive Control ⁵	LOT 62641156 Negative Control ⁶
CTX-M		Detected ⁷	Not Detected ⁸
IMP		Detected ⁷	Not Detected ⁸
KPC		Detected ⁷	Not Detected ⁸
NDM		Detected ⁷	Not Detected ⁸
OXA		Detected ⁷	Not Detected ⁸
VIM		Detected ⁷	Not Detected ⁸

- streck.com/IFU
- 1

Quality Control for the Luminex VERIGENE[®] BC-GN Test
Contrôle de la qualité du Luminex VERIGENE BC-GN Test / Qualitätskontrolle für das Luminex VERIGENE BC-GN Test / Controllo di qualità per il Luminex VERIGENE BC-GN Test / Control de calidad de Luminex VERIGENE BC-GN Test / Kvalitetskontroll för Luminex VERIGENE BC-GN Test
- 2

The product may be used until expiration date.
Une fois ouvert, le produit peut être utilisé jusqu'à la date de date de péremption. / Einmal geöffnet kann das Produkt bis zum angegebenen Ablaufdatum verwendet werden. / Una volta aperto, il prodotto può essere utilizzato fino alla data di scadenza. / Una vez abierto, puede utilizarse el producto hasta la fecha de caducidad. / När produkten har öppnats kan den användas fram till utgångsdatumet.
- 3

Gram-Negative Bacteria
Bactérie Gram-négative/ Gramnegative Bakterien/ Batteri Gram-negativi/ Bacterias Gramnegativas/ Gramnegativa bakterier
- 4

Target
Cible/ Ziel / Target/ Objetivo/ Mål
- 5

Positive Control
Contrôle positif/ Positive Kontrolle/ Controllo positivo/ Control Positivo/ Positiv kontroll
- 6

Negative Control
Contrôle négatif/ Negative Kontrolle/ Controllo negativo/ Control Negativo / Negativ kontroll
- 7

Detected
Détecté / Entdeckt / Rilevato / Detectado / Detekterad
- 8

Not Detected
Non détecté / Nicht entdeckt / Non rilevato / No detectado / Ej detekterad
- 9

Antimicrobial Resistance Genes
Gènes de résistance aux antimicrobiens / Antimikrobielle Resistenzgene / Geni di resistenza antimicrobica / Genes con resistencia antimicrobiana / Antimikrobiella resistensgener
- 10

Gene
Gène/ Gene/ Gene/ Gen/ Gen

All product names, logos, brands, and marks are property of their respective holders.
Tous les noms, logos, marques et labels de produits sont la propriété de leurs propriétaires respectifs. / Alle Produktnamen, Logos, Marken und Zeichen sind Eigentum ihrer jeweiligen Besitzer. / Tutti i nomi dei prodotti, i loghi, i marchi e le marche sono di proprietà dei rispettivi titolari. / Todos los nombres de productos, logotipos, marcas comerciales y otras marcas son propiedad de sus respectivos propietarios. / Alla produktnamn, logotyper, varumärken och märken tillhör respektive innehavare.