



June 2017

Streck is pleased to inform you of a change in the storage temperatures for empty Cell-Free DNA BCT®. This change has been internally validated and can be implemented with the lot numbers listed below.

New Storage Guidelines for Empty Tubes

When stored at 2 °C - 30 °C, empty Cell-Free DNA BCT is stable through expiration date.

Short-term storage from 2 °C - 40 °C is acceptable for empty Cell-Free DNA BCT for up to 14 days.

Description	Catalog Numbers	Lot Number	Expiration Date	
Cell-Free DNA BCT		5210	2017-07-29	
		5308	2017-11-04	
		6027	2018-01-27	
		218961	6111	2018-04-20
		218962	6115	2018-04-20
		218992	6195	2018-07-13
		230242	4255	2016-08-20
		230243	5054	2017-01-14
		6279	2018-10-05	
		6363	2018-12-28	
7079	2019-03-20			
Cell-Free DNA BCT CE	218996	5252	2017-09-09	
		6083	2018-03-23	
		6251	2018-09-07	
		7025	2019-01-25	
		230244		
230250				
230251				

If you have additional questions, please contact Streck Technical Services at 800.843.0912.

Best Regards,

Justyna Tesmer  
Product Manager