MDx-Chex® for MVP

streck.com/IFU

STRECK

Quality Control for MVP Assays

350861
2026-03

For Research Use Only.
Not for use in diagnostic procedures.

2027-11-08

The product may be used until expiration date.

Expected Results

Bacteria		
Target	LOT **63131201 Positive Control**	LOT **63131211 Negative Control**
Bacterial Vaginosis	Positive	Negative
Yeast		
Candida Group	Detected	Not Detected
Candida-glab-krus	Detected	Not Detected
Parasites		
Trichomonas vaginalis	Detected	Not Detected

All product names, logos, brands, and marks are property of their respective holders.

Streck
7002 S. 109 Street
La Vista, NE 68128 USA