MDx-Chex® for BCY

IVD

STRECK

Quality Control and Verification Panel for Diasorin LIAISON PLEX® Yeast Blood Culture Assay

350829
2026-05

streck.com/IFU

2027-11-04

Rx Only

The product may be used until expiration date.

Expected Results

Yeasts		
Target	LOT **63341170 Positive Control**	LOT **63341181 Negative Control**
Candida albicans	Detected	Not Detected
Candida dubliniensis	Detected	Not Detected
Candida glabrata	Detected	Not Detected
Candida krusei	Detected	Not Detected
Candida parapsilosis	Detected	Not Detected
Candida tropicalis	Detected	Not Detected
Candida auris	Detected	Not Detected
Candida kefyr	Detected	Not Detected
Candida lipolytica	Detected	Not Detected
Candida lusitaniae	Detected	Not Detected
Candida guilliermondii	Detected	Not Detected
Cryptococcus neoformans/ gattii	Detected	Not Detected
Candida famata	Detected	Not Detected
Candida haemulonii/ duobushaemulonii	Detected	Not Detected

All product names, logos, brands, and marks are property of their respective holders.

Streck
7002 S. 109 Street
La Vista, NE 68128 USA